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Please store this document near the device for later use!

MedizinTechnik
... for a better life

English

Operating Instructions

Atmolit 26

320.0200.B
ON.09

	Page
1.0 Introduction	3 - 5
1.1 Notes on the operating instructions	3
1.2 Function	4
1.3 Extents of supply.....	4
1.4 Transport and storage.....	4 - 5
1.5 Explanation of symbols.....	5
2.0 Safety advice	6
3.0 Intended use	7
4.0 Setting up and starting up	8
5.0 Operation	9
6.0 Notes on cleaning and care	10 - 13
6.1 Basic information	10
6.1.1 Cleaning the unit's surfaces.....	11
6.2 Evacuating the collection jar	11
6.2.1 Cleaning the collection jar parts.....	11
6.3 Recommended disinfectants for instruments	12
6.4 Recommended disinfectants for surfaces.....	12
6.5 Application of the disposable	13
bacterial filter plate	
6.6 Reprocessing	13
7.0 Maintenance and servicing	14
7.1 Basic information	14
7.2 Change of the disposable	14
bacterial filter plate	
8.0 Trouble-shooting.....	15
9.0 Accessories, consumables,.....	16 - 17
spare parts	
9.1 Accessories.....	16
9.2 Consumables	17
9.3 Spare parts	17
10.0 Technical specifications.....	18-19
11.0 Disposal	20
Declaration of conformity	

1.1 Notes on the operating instructions

- These operating instructions apply for following devices:
 - Atmolit 26 / G.....REF 320.0200.0
with glass jar 1 l
 - Atmolit 26 / K.....REF 320.0201.0
with TPX plastic jar 1.25 l
 - Atmolit 26 / R.....REF 320.0202.0
with disposable Receptal®jar 1.5 l
 - Atmolit 26 / G-2REF 320.0203.0
with glass jar 2 l
 - Atmolit 26 / G-SREF 999.0200.0
with glass jar 1 l and
additional rinsing jar

In combination with the accessories mentioned in chapter 9.0 as well as for all identically constructed devices with special voltages.

Further information, accessories, consumables and spare parts are available from:

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Fax:
+ 49 7653 689-392 (Sales Germany)
+ 49 7653 689-391 (Sales Worldwide)

atmos@atmosmed.de
<http://www.atmosmed.de>

- These operating instructions contain important notes on how to operate the Atmolit 26 safely, correctly and effectively. Therefore, they are intended not only for new operating personnel to be instructed in its use, but also for use as a reference manual. Following the operating instructions will help you to avoid risks, reduce repair costs and down-time and increase reliability and service-life of the equipment. **These operating instructions must always be kept available near the equipment.** Prior to first use please read chapter 2.0 "For your safety", to understand and avoid possible dangerous situations.

The basic principles are:

Judicious and careful work provides best protection against accidents!

Operational safety and readiness for use depend not only on your capabilities, but also on **care and maintenance** given to the Atmolit 26. For this reason regular cleaning and service work are a must. Major maintenance and repair work may be carried out only by expert personnel authorized by ATMOS. In case of repairs you should insist that original spare parts only are used to continue the warranty and ensure the operational safety and readiness of the unit.

- The product Atmolit 26 bears CE-0124 marking according to the EEC guideline of the council for medical products 93/42/EEC and meets the basic requirements of annex I of this guideline.
- The quality management system used at ATMOS has been certified for quality management in accordance with the relevant international standards.
- For authorised servicing, ATMOS provides servicing instructions with detailed descriptions of circuits, setting instructions and servicing information.
- These operating instructions correspond to the design of the device and the status of basic safety engineering standards on going to press.
- Reproduction of these instructions – even in part – only with the written permission of ATMOS.

Abbreviations/symbols:

- Identifies a list
 - Subdivision of a list/activity.

The recommended sequence must be adhered to in each case!

 Important advice!

 Description of the effect of an activity.


Fig. 1. Atmolit 26

- ① Suction hose
- ② Attachment (for suction hose)
- ③ Filter attachment
- ④ Silicone connecting hose
- ⑤ Vacuumgauge
- ⑥ Vacuum adjustment
- ⑦ Hose attachment (on the unit)
- ⑧ Collection jar
- ⑨ Green signal on filter attachment (③) and attachment on the unit (⑦)

1.2 Function

The Atmolit 26 is a very handy small suction unit. It is driven by an electromotive, maintenance-free diaphragm-type pump. During operation, the pump generates a vacuum within the hose system and the collection jar, thus sucking off secretions or fluids (e.g. by means of a suction catheter). The fluid is gathered in the collection jar. A mechanical overflow safety (on the inner part of the collection jar lid) avoids penetration of secretion into the pump head. The final vacuum and, following, the air-flow rate can be adjusted by means of the fine control and the vacuumgauge. The unit is equipped with a rechargeable battery (accumulator). Integrated microprocessor-based technology assures safe charging of the battery; overcharging is thus impossible. An overtemperature stop controlled by electronics avoids overheating of the unit. A disposable bacterial filter plate integrated in the lid of the collection jar prevents bacteria and liquid from penetrating into the pump.

1.3 Extents of supply

- Prior to dispatch, this ATMOS device was subjected to an extensive functional test and has been carefully packed. Nevertheless, please compare the contents of the shipment on completeness immediately upon receipt (see delivery note).

1.4 Transport and storage

- The transport of the device may be effected only in a dispatch carton upholstered and offering sufficient protection.
- Please document and report damages in transit immediately. For complaints or return deliveries, please use the enclosed form **QD 434**.

1.4 Transport and storage (contin.)

- Ambient conditions:

Transport / storage: -30...+50°C;
 20...80 % humidity
 without condensing
 air pressure 700...1060 hPa

Operation: +5...+35°C;
 20...80 % humidity
 without condensing
 air pressure 700...1060 hPa

- The unit must be allowed to stand for up to six hours at room temperature prior to starting up for the first time following transport at temperatures below freezing. The unit may **not** be operated if it has **not** acclimatised as this might damage its diaphragms.

1.5 Explanation of symbols

Important advice !

Symbol
"Caution: pay attention to operating instructions"

This symbol shall draw your attention to the safety information contained in the operating instructions. It serves as a reference to particular notes which are important for the application of Atmolit 26.

Protection class II

Unit fuse

Insulating class type B (body)

The CE sign shows that this product meets the appropriate requirements of the EU guidelines.

SN

Serial number

REF

Ordering number

Date of production

- Prior to starting up the Atmolit 26, read these operating instructions carefully.
- ATMOS cannot guarantee perfect functioning neither will it be liable for damage to people or property if:
 - Any non-original ATMOS parts are used,
 - the user instructions given in this manual are not followed exactly or are disregarded,
 - Assembly, resetting, alterations, extensions and repairs are not carried out by people authorised by ATMOS.
- No warranty rights shall exist in the event of damage or failure caused by the use of non-ATMOS accessories or non-ATMOS consumables.
- The safety standard of the Atmolit 26 corresponds with recognized medical technical regulations and the directions of the law relating to medical products.
- Only persons instructed in medical use may apply the Atmolit 26 to patients.
- The unit may not be started:
 - If cables or plugs are defective,
 - if it has been dropped down before,
 - if obvious defects might restrict safe operation.Prior to returning the device for repair, clean it.
- Carry out a visual inspection of the unit prior to each use including hoses, collection jar and connection cable.
- The unit may not be operated in splash water range and in locations where there is a danger of explosion (zones M and G).
- Pay attention to the ambient conditions described in chapter 1.4 Transport and storage.
- Never plunge the unit into water, not even when it is switched off.
- Do not allow any liquid to get into the unit. If liquid has penetrated the unit, it may not be operated again until it has been checked by the customer service centre.
- Inserting metallic subjects through the louvres in the interior of the unit might involve danger to life.
- The Atmolit 26 has been designed for aspirating body fluids in medical ranges. Never remove explosive gases and inflammable or corrosive fluids.
- The unit must be set up on a firm, level surface. The switched-on unit might get overheated if it is placed on an uneven surface (e.g. mattress, cushion, padded seat etc.).
- The main voltage specified on the type plate must match the power supply system.
- Never connect the unit to defective power sockets or extension cables.
Avoid moisture on plug and switches.
- Use transparent hoses exclusively.
- This suction unit may not be applied without disposable bacterial filter plate.
- The unit, collection jar, mains cable, accessories, connection cables and hoses must be checked for damage prior to starting up. Damaged cables and hoses must be replaced immediately. **Prior to use, check the unit functions.**
- The suction hose must never come into direct contact with the application site. A suction catheter, attachment or a medical aspiration set must always be connected to the hose.
- Sterile packed parts may no longer be used if their packing was damaged during transport or storage ⇨
Danger of infection for the patient.
- Always remove the plug from the wall socket first in order to disconnect the unit from the mains. Only then may the connecting cable be disconnected from the unit:
 - Before cleaning the unit,
 - before the collection jar is evacuated,
 - before leaving the room.Never pull at the cable !
Never touch the plug or cable with wet hands.

- The Atmolit 26 is used in medical ranges for the suction of secretions and body fluids.
- **Special medical ranges for applying this unit are:**
 - Ambulant and home care of old people for spontaneous suction of mucus from the respiratory tract;
 - first aid posts, ambulance stations, wards, old people's and nursing homes;
 - military infirmaries, field hospitals;
 - persons who require on-going removal of secretion from the Trachea.
- The suction hose must never come into direct contact with the application site. A suction catheter, attachment or a medical aspiration set must always be connected to the hose.
- The Atmolit 26 may not be used:
 - for drainages in low vacuum range (e.g. thoracic drainages) and for suction procedures outside medical ranges;
 - for the suction of inflammable, corrosive or explosive fluids/gases.


Fig. 2. Rear view of the Atmolit 26

- ① Power attachment
- ② Power switch
- ③ Fuse support

- The Atmolit 26 is delivered ready for use.
- Lift the unit out of the cardboard.
Check whether the voltage values on the data plate correspond with the inbuilding voltage.
- Prior to first operation, pay attention to the safety information in chapter 2.0.
- The unit must be allowed to stand for up to six hours at room temperature prior to starting up for the first time following transport at temperatures below freezing. The unit may **not** be operated if it has **not** acclimatised as this might damage its diaphragms.
- Join suction hose (①, fig. 1, page 4) with the attachment (②, fig. 1, page 4).
- Connect silicone hose (④, fig. 1, page 4) to filter attachment (⑤, fig. 1, page 4) and attachment on the unit (⑦, fig. 1, page 4) (pay attention to green signal ring).
- Connect the power cable (①, fig. 2) to the power attachment and the power plug to the local power socket.
- Start the unit with the switch (②, fig. 2).

Changing fuse (fig. 2)

Loosen the fuse support by compressing the lower spring clip and pull it out (use a 2 mm flat screwdriver, if necessary). Replace the fuse (see fuse values in chapter 10.0 Technical specifications) and insert the fuse support again.


Fig. 3. Vacuumgauge


Fig. 4.

- ① Suction hose
- ② Finger tip
- ③ Suction catheter

AUXILIARY AIR VENT OPEN = suction procedure is inter-rupted (e.g. when leading in the catheter)

AUXILIARY AIR VENT CLOSED WITH THE FINGER = suction

i This suction unit may not be applied without disposable bacterial filter plate as, otherwise, the unit might get contaminated or the pump might get destroyed.

i **Attention: Suction procedures in the respiratory tract may only be implemented after appropriate instruction by hospital or special staff.**

If the Atmolit 26 is used e.g. for suction of mucus from the upper respiratory tract, proceed as follows:

- Adjust your desired maximum vacuum by closing the suction hose opening with the finger (①, fig. 4) the vacuum is then generated. Open the regulating valve (⑥, fig. 1, page 4) until the vacuumgauge shows the desired vacuum value (fig. 3).
- Choose a suction catheter of the right size (e.g. Unoplast catheters ③, fig. 4, which are available from ATMOS in 3 different sizes) or a suction instrument which is only available from specialized dealers.
- Join the suction hose (①, fig. 4) and the suction catheter (③, fig. 4) by means of the finger tip (②, fig. 4).
- Then, lead in the catheter in the same way as shown by your hospital staff and start the suction procedure. Never start suction procedures without having been instructed by hospital staff.
- Control the suction procedure with the finger tip.
- Rinse suction catheter and suction tube with clean water after every suction procedure. We recommend to use the rinsing bottle (REF 000.0504.0) with lid (REF 000.0504.1). Thus, you can carry clean water along with you.

i **Make sure that the collection jar is evacuated in time. As soon as the jar is half-filled, it must be emptied (this principle proves right in all application ranges).**

- When the maximum level is exceeded, the overflow safety reacts and suction is stopped. Empty the jar as described in section 6.2.

i **If secretion has penetrated the pump, the Atmolit 26 has to be maintained by a service technician.**

6.1 Basic information

- Remove the mains plug from the socket before starting cleaning and disinfection!
- The measures described for cleaning and disinfection or sterilisation are not intended to replace the regulations applicable for operation in each case!
- Always comply on principle with the relevant manufacturer's information regarding concentration and notes for use!
- **Attention:** The lid parts and silicone hoses might get dyed by some disinfectants; a fact which does not take effect on the attributes of the materials. These silicone and collection jar parts may also be boiled (for more than 10 min).
- Do not use any
 - **Disinfectants containing organic or inorganic acids or bases since these may cause corrosion damage.**
 - **Disinfectants containing chloramines, phenol derivatives or anionic tensides since these may cause stress cracks in the plastics used.**
- Basically, all parts which come into contact with secretion must be cleaned, disinfected or sterilized after every use.
- Silicone hoses are autoclave safe (up to 136°C). These parts can also be immersed in commercial disinfectants (see page 12).
- The disassembled components of the secretion container (container, lid, overflow protection) can be autoclaved (up to 136°C). It is also possible to immerse these parts into standard disinfectant solutions (see page 12).

6.1.1 Cleaning the unit's surfaces

If liquid has got into the unit, it may not be operated again until it has been checked by the authorised customer service centre.

- The surfaces of the Atmolit 26 are resistant to all the surface disinfectants listed on page 12.
- The unit itself can be wiped off with a moist cloth (not wet).

6.2 Evacuating the collection jar

- Remove connecting hoses (suction hose and connecting hose on the unit, fig. 5).
- Lift out collection jar with lid.
- Remove collection jar lid.
- The collection jar can then be evacuated.
- Properly dispose of the sucked material.


Fig. 5.

6.2.1 Cleaning the collection jar parts

- Collection jar and collection jar lid may be rinsed under running water or cleaned in an automatic cleaning device.

Using the cleaning agent Neodisher AN (manufactured by Dr. Weigert, Hamburg) cleaning in a special dish-washer is also possible. The disposable bacterial filter plate (fig. 6, page 13) must be removed before cleaning.

6.3 Recommended disinfectants for instruments

Disinfectant	Contents	(in 100 g)	Manufacturer
GIGASEPT FF (Application concentrate)	Succindialdehyde Dimethoxytetrahydrofurane Corrosion inhibitors Non-ionic tensides and fragrances	11.0 g 3.0 g	Schülke & Mayr, Norderstedt
Sekusept PLUS ¹ (Application concentrate)	Glucoprotamine Non-ionic tensides Solvents, complexing agents	25.0 g	Ecolab, Düsseldorf
Mucozit-T (Application concentrate)	Bis(3-aminopropyl)laurylamine Alkyl dimethylbenzyl ammonium chloride Cocosporylendiamin-1,5-guanidiniumacetat	8.0 % 19.0 % 7.0 %	Merz & Co., Frankfurt/Main

6.4 Recommended disinfectants for surfaces

Disinfectant	Contents	(in 100 g)	Manufacturer
TERRALIN (Application concentrate)	Benzalkonium chloride Phenoxypropanoles	20.0 g 35.0 g	Schülke & Mayr, Norderstedt
QUATOHEX (Application concentrate)	Didecyl dimethyl ammonium chloride Benzalkonium chloride Bi-guanidinium acetate Polymer biguanide Active cleaning substances	14.0 g 10.0 g 7.5 g 0.5 g	Braun, Melsungen
Incidin Plus (Application concentrate)	Glucoprotamine Non-ionic tensides Solvents, complexing agents	26.0 g	Ecolab, Düsseldorf
Pursept-A (Disinfectant spray or disinfectant cloths)	Ethanol Glyoxal QAV	38.9 g 0.1 g 0.05 g	Merz & Co., Frankfurt/M.

Discolouration may result if disinfectants containing aldehydes and amines are used on the same object.

6.5 Application of the disposable bacterial filter plate

i This suction unit may not be applied without bacterial filter !
The bacterial filter plate **cannot be cleaned or reused again !**

- The collection jar lid contains a disposable bacterial filter plate (2, fig. 6), which protects the interior of the unit against contamination and stops fluids due to its hydrophobic attribute. Together with the mechanical overflow safety, this filter avoids the sucking-in of secretion into the interior and, thus, the failure of the unit. The basic element of the filter material is a porous teflon layer which is put on a supporting material (polyester fleece). The pore size of the teflon material is 1µm.
- The bacterial filter plate is a disposable part, i.e. it is used until it is worn resp. until, for hygienic reasons, it can no longer be used. The disposable bacterial filter plate is worn when the vacuumgauge on the suction unit shows a difference vacuum value between -30 kPa and -40 kPa while the system is in closed and the suction hose in open condition. Hygiene is no longer guaranteed when the disposable bacterial filter plate is wet as it might then be a medium for germs and might affect functioning of the unit. The disposable bacterial filter plate must also be changed when it is apparently contaminated.
- If the suction unit is used in one and the same patient the disposable bacterial filter plate must be replaced every other day at the latest. If the suction unit is used in several patients the bacterial filter plate must be replaced daily.


Fig. 6. Disposable bacterial filter plate

- 1 Collection jar lid
- 2 Disposable bacterial filter plate
- 3 Bacterial filter cover

Disposable bacterial filter plates,REF 320.0065.0
25 pcs.

Disposable bacterial filter plates,REF 320.0066.0
100 pcs.

6.6 Reprocessing


The way the suction device is used determines its reliability and safety. These hygiene measures described in the last chapter are indispensable for protecting the patient and the user and for maintaining a safe and reliable suction device.

These measures do not replace a reprocessing, performed by the manufacturer or by any certified ATMOS partner before re-using the device in a new patient. Before using the suction device with a new patient it has to be reprocessed by ATMOS or a certified ATMOS partner in order to protect the user/patient. How can one realise that the suction device is contaminate?

- 1.) Perform a visual inspection of the connecting hose and of the connecting adapter (at the device, see fig. 1). In the case one or both parts are soiled or damp, the device is oversucked and therefore contaminated.
- 2.) In the case the cover of the bacterial filter (fig. 2) is damp or soiled, the device is oversucked and therefore contaminated.
- 3.) If the insulating mat or the sound absorber are wet or soiled, the device is oversucked and therefore contaminated. For this inspection you have to remove the cover from the sound absorber (fig. 3).

In case you realise one of the aforesaid incidents, it is necessary to have the device repaired by ATMOS or by a certified ATMOS Partner.

If any reservations exist regarding the hygiene condition of the device, please send the device to ATMOS or a certified partner for inspection.

7.1 Basic information

- Carry out a visual inspection of the unit prior to each use including hoses, collection jar and connection cable. **Damaged cables and hoses must be replaced immediately.**
- Maintenance or opening and repair of the Atmolit 26 (with the exception of the cleaning work described in these operating instructions) may only be carried out by ATMOS or a specialist authorised by ATMOS. In this case, attention should be paid to the protective technical and hygiene measures, the notes on safety plus the descriptions in the servicing instructions for the Atmolit 26.
- For repair, this device can be returned to ATMOS.
- Before returning the device for repair, clean and disinfect it.
- ATMOS cannot guarantee perfect functioning neither will it be liable for damage to people or property if:
 - Any non-original ATMOS parts are used,
 - the user instructions given in this manual are not followed exactly or are disregarded,
 - assembly, resetting, alterations, extensions and repairs are not carried out by people authorised by ATMOS.
- No warranty rights shall exist in the event of damage or failure caused by the use of non-ATMOS accessories or non-ATMOS consumables.
- Pay attention to regulations and instructions valid for the respective application range.

7.2 Change of the disposable bacterial filter plate

- Remove silicone hoses on collection jar lid.
- Open bacterial filter cover by turning it to the left. Simultaneously hold silicone ring and attachment of collection jar lid.
- Lift out the disposable bacterial filter plate and replace it by a new one.

The smooth side of the disposable bacterial filter plate must always show in direction of the collection jar. The blue mark must show in direction of the bacterial filter lid.

Please make sure that one end of the silicone hose is connected to the attachment on the unit and the other end to the filter attachment (③, fig. 1, page 4) (green rings).

Attention: Failure of the unit will occur if the silicone hose is connected to the attachment for suction hose (②, fig. 1, page 4).

Prior to dispatch, the Atmolit 26 was subjected to an extensive functional test. If, nevertheless, a failure should appear, you may possibly clear it yourself if you follow these notes:

Problem	Possible causes	Remedy
<ul style="list-style-type: none"> ● Unit does not start 	<ul style="list-style-type: none"> – Loose power plug – No power voltage – Defective fuse 	<ul style="list-style-type: none"> – Check connection to supply socket – Check inbuilding fuse – Check whether power plug on the unit is loose – Replace fuse
<ul style="list-style-type: none"> ● Insufficient performance 	<ul style="list-style-type: none"> – Leakages within the hose system or in collection jar lid – Clogged filter (vacuumgauge indicates vacuum) 	<ul style="list-style-type: none"> – Check collection jar lid and hose system, replace sealing ring on collection jar lid, if necessary – Replace filter
<ul style="list-style-type: none"> ● No suction 	<ul style="list-style-type: none"> – Float of overflow safety closes the collection jar inlet – Blocked filter – Suction hose is connected to filter cover – Secretion or blood has been sucked in and valve plates of the pump are contaminated 	<ul style="list-style-type: none"> – Check filling level of collection jar; evacuate jar, if necessary; clean overflow safety and check free movement of the float – Replace filter in collection jar lid and check correct hose connection – Change silicone hoses on collection jar lid (green rings) – Unit has to be returned for repair


Fig. 7. Atmolit 26
 ① Suction hose
 ② Connecting hose
 ③ Top to the handle
 ④ Collection jar (glass or TPX plastic)
 ⑤ Collection jar lid, complete


Fig. 8. Collection jar lid
 ① Bacterial filter cover
 ② Sealing ring for bacterial filter
 ③ Disposable bacterial filter plate with blue marking
 ④ Insert of collection jar lid
 ⑤ Gasket of collection jar
 ⑥ Sleeve of overflow safety
 ⑦ Float, closed side must show upwards


Fig. 9.
 ① Suction hose
 ⑥ Finger tip
 ⑦ Suction catheter

9.1 Accessories

Fig.	Description	REF
④, Fig. 7	grad. collection jar 1 l, glass	000.0057.0
No Fig.	grad. collection jar 2 l, glass	000.0618.0
④, Fig. 7	grad. collection jar 1.25 l, (TPX) vapour-sterilizable up to 136°C	000.0544.0
No Fig.	Retaining bracket for Receptal® container	320.0020.0
No Fig.	Receptal® container for suction bags	310.0221.0
No Fig.	Receptal® suction bags 1.5 l with integr. overflow valve filter 50 pcs.	310.0222.2
No Fig.	Medi-Vac® support	320.0140.0
No Fig.	Wall- / rail unit fixture	320.0120.0
No Fig.	Trolley, on 4 castors self-assemblage	320.0070.0
Attention: Trolley and small suction unit cannot be screwed up in standard version! If required, we will deliver a special design of the suction unit.		
No Fig.	Basket with standard rail holder	320.0075.0
Unoplast suction catheters "Optimal", straight, central opening, 2 small lateral openings, length: 50 cm, each packed separately and sterile in packing lots of 100 pcs.		
⑦, Fig. 9	Size: charrière 12	000.0294.0
⑦, Fig. 9	charrière 14	000.0295.0
⑦, Fig. 9	charrière 16	000.0296.0
i	Sterile packed parts may no longer be used if their packing was damaged during transport or storage ⇨ Danger of infection for the patient.	
No Fig.	Operating instructions	320.0200.A

9.2 Consumables

Fig.	Description	REF
⑤, Fig. 7	Collection jar lid, complete.....	320.0060.0
④, Fig. 8 Insert of collection jar lid	320.0011.0
⑤, Fig. 8 Gasket for insert of collection jar lid	320.0013.0
①, Fig. 8 Bacterial filter cover	320.0012.0
⑥, Fig. 8 Sleeve of overflow safety.....	320.0010.0
⑦, Fig. 8 Float.....	320.0015.0
No Fig. Rinsing jar 250 ml	000.0504.0
No Fig. Lid for rinsing jar	000.0504.1
No Fig. Medi-Vac® external container 1 l	312.0473.0
No Fig. Medi-Vac® internal container 1 l. 50 pcs.	312.0474.0
③, Fig. 8 Disposable bacterial filter plates, 25 pcs.	320.0065.0
③, Fig. 8 Disposable bacterial filter plates, 100 pcs.	320.0066.0
②, Fig. 8 Gasket for filter cover.....	320.0016.0
No Fig. Suction hose, silicone, Ø 6 mm, 1.30 m, (136°C)	000.0013.0
No Fig. Suction hose, Ø 6 mm, vapour-sterilizable up to 136°C, price / m	006.0009.0
⑥, Fig. 9 Hose connector (finger tip), sterile, 1 pcs.	000.0347.0
⑥, Fig. 9 Hose connector (finger tip), sterile, 100 pcs.	000.0347.1
No Fig. Suction hose, disposable, Ø 6 mm, 0.25 m, sterile, 10 pcs.	006.0058.0
No Fig. Suction hose, disposable, Ø 6 mm, 1.30 m, sterile, 10 pcs.	006.0057.0
No Fig. Suction hose, disposable, Ø 6 mm, 2.10 m, sterile, 50 pcs.	006.0059.0

Sterile packed parts may no longer be used if their packing was damaged during transport or storage
 ⇨ Danger of infection for the patient.

No Fig. Power cable, 2-wire	008.0647.0
---------	---------------------------------	------------

9.3 Spare parts

Fig.	Description	REF
③, Fig. 2 Fuse T 0,500 A/H	008.0719.0
②, Fig. 7 Connection hose, silicone, Ø 6mm, 28 cm	320.0049.0
③, Fig. 7 Top of the handle	320.0005.0
No Fig. Service valve	320.0130.0

Air flow rate of pump	26 ± 2 l/min
Max. vacuum at sea level	-80 kPa*
Vacuum readout	-1...0 bar (± 25 mbar) (mm Hg; kPa)*
Additional air regulation	mechanical regulating valve
Collection jar	1 l / 2 l jars made of glass or 1.25 l plastic (TPX) or 1.5 l Receptal [®] -system
Suction hose	∅ 6 mm, 1.30 m length
Voltage	230 V~ 50/60 Hz Special voltage: 115 V~ 50/60 Hz
Current input (max.)	approx. 0.3 A for 230 V~ approx. 0.6 A for 115 V~
Power consumption	approx. 65 W
Power cable	2 m
Fuse	slow-blow 500 mA/H for 230 V~ slow-blow 1.2 A/H for 115 V~
Operating time	Interrupted use over approx. 60 min., cooling period approx. 60 min., depending on ambient temperature
Protective earth conductor resistance	—
Earth leakage current	—
Enclosure leakage current	N.C. < 0.1 mA
Patient leakage current	—
Heat emission	65 J/s
Noise level	55.4 dB (A) @ 1m (acc. to ISO 7779)
Ambient conditions	-30...+50°C
Transport/storage	20...80 % humidity, non-condensing air pressure 700...1060 hPa
Operation	+5...+35°C 20...80 % humidity, non-condensing air pressure 700...1060 hPa
Dimensions HxBxT	206 x 362 x 174 mm
Weight	4 kg
Protection class (EN 60601-1)	II
Degree of protection	Type B 
Protection category	IP00
Classification acc. to Annex IX EEC directions 93/42/EEC	Ila
CE marking	CE 0124
Rules applied	EN 60601-1: 1990 +A1: 1993 + A2: 1995 EN 60601-1-2: 1993 (EMV / EMC) EN 30993: 1994
UMDNS-Code	10-219
REF	320.0200.0 (with 1 l jar made of glass) 320.0203.0 (with 2 l jar made of glass) 320.0201.0 (with 1,25 l TPX-jar) 320.0202.0 (with 1,5 l Receptal [®] -system)

*1 bar ≙ 750,06 mm Hg ≙ 1000 hPa / depends on daily atmospheric pressure

Canadian Classification

Device Group

PNC

Risk Class

Description

Ear, Nose and Throat

77QBW

2

ASPIRATOR, TRACHEAL

Issue of Technical Specifications: 01.04.2005

- The Atmolit 26 is not comprised of any hazardous materials.
- The materials of the housing can be recycled completely.
- Prior to disposal, device and accessories must be decontaminated.
- The materials are to be separated carefully.
- Pay attention to country-specific regulations for disposal (e. g. waste incineration).

Cleaning and servicing plan for ATMOS LC 16, Atmolit, ATMOPORT

start date::

name of the item:

serial number:


--	--	--

day	cleaning secretion container	cleaning container lid	cleaning of the housing	exchange of bacterial filter	exchange of fingertip	exchange of suction hose, 1.3 m	exchange of connecting hose	cleaning/exchange performed	
								name	signature
				exchange*	exchange*	exchange*	exchange*		
					exchange				
				exchange	exchange				
					exchange				
				exchange	exchange				
					exchange				
				exchange	exchange				
					exchange				
				exchange	exchange				
					exchange				
				exchange	exchange				
					exchange				

daily, respectively after each use

daily, respectively after each use

daily, respectively after each use

* Before first time operation of a brand new device, respectively a reprocessed device, an exchange is unnecessary!

Special notes:

Before operating the suction device, the user has to make sure that the device functions and is in good order and condition. The user has to observe the instructions in the operating manual as well as all other safety-related and maintenance information enclosed.

For cleaning and disinfection only agents which are recommended by the manufacturer, may be used.

Only sterile, single-use suction catheters may be used for suctioning. They have to be exchanged before each suction process. During use utmost attention to hygiene (e.g. disinfection of hands, wearing single-use gloves) is indispensable. After each use the secretion container and the hose must be rinsed thoroughly with water. During storage the contamination of the device

and other products must be avoided. The intervals stated in the list are non-binding guide values. Depending on the use shorter intervals may be necessary.

For each patient a new or a reprocessed suction device must be used. Otherwise there is high and acute danger of infection for the patient, the user and any third person!

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Important notes

General information


The way the suction device is used determines its reliability and safety. These hygiene measures are indispensable for protecting the patient and the user and for maintaining a safe and reliable suction device.

These measures do not replace a reprocessing, performed by the manufacturer or by any certified ATMOS partner before re-using the device on a new patient.


This cleaning and servicing plan as well as the relevant notes result from many years of experience. Depending on the use and the user's experience shorter intervals may be necessary.

ATMOS recommends the following sets of consumables:


320.0125.0 für LC 16, Atmolit 16N, Atmolit 26

312.0615.5 für ATMOPORT, ATMOPORT S, ATMOPORT N

suction catheter, length: 50 cm, 100 pieces

- ☞ Disconnect the mains plug from the socket before commencing with cleaning and disinfection!
- ☞ Please observe the notes in the operating instructions, especially regarding the recommended agents.
- ☞ All parts (except bacterial filter and device) are autoclave-safe up to 130°C.

Cleaning of the secretion container

Please empty the secretion container after each use, rinse it thoroughly with warm water and clean it with washing-up liquid. Tenacious

contaminations can be removed with a standard bottle brush.

Cleaning of the container lid

The bacterial filter must be removed before cleaning, please use single-use gloves or tweezers. Please demount the container lid after each use and rinse it thoroughly. The lid must be absolutely dry before reuse. Please pay attention to a correct function of the overflow safety when mounting the lid.

Bacterial filter

The bacterial filter prevents penetration of micro organisms and secretion into the device, respectively blowing out from it and is therefore a protection for the user and the device. For hygienic reasons a weekly exchange is recommended. If the maximum vacuum is adjusted, the suction hose is open and the vacuum gauge shows a basic vacuum of > -0.3 bar, then the bacterial filter must be exchanged immediately. In case of contamination the filter must also be replaced. In order to increase the service life of bacterial filters, it is recommended to empty the secretion container when it is half-full. Always use the original ATMOS bacterial filter.

☞ Never use the suction device without bacterial filter!

Hose connection/fingertip

The fingertip connects the suction hose to the suction catheter. The fingertip is in continuous contact with secretion and is difficult to clean, therefore we recommend an exchange every 2 or 3 days.

Suction Hose

The suction hose conducts the secretion from the suction catheter to the container. In order to prevent secretion from drying, the hose must be thoroughly rinsed with clear water after each use. The water can be sucked into the secretion container. Please fill the secretion container only half. Frequent cleaning and disinfection/sterilisation may discolour and embrittle the hose. Therefore, a monthly exchange of the suction hose is recommended.

Connecting hose

This hose is the connection between secretion container and the vacuum pump. For hygienic reasons a monthly exchange of the hose is recommended.

Coding rings

On the LC 16, Atmolit 26 and Atmolit 16 N, the green coding rings mark the positions where the connecting hose is connected to the container lid respectively bacterial filter housing and the nipple at the device.

Cleaning of the device (housing)

When the device is contaminated but at least once per week the housing must be wiped off with a moist (but not wet) cloth. A weekly disinfection is recommended.

☞ Never irrigate the device with water and never emerge it into any liquid.

Cleaning/disinfection

To improve the cleaning effect, standard washing-up liquid can be added to the warm water. In the case of tenacious contamination the parts should be steeped in water for a length of time or they may be removed with a soft brush or cloth. After thorough cleaning, container, fingertip and hoses can be disinfected with a disinfection agent (see operating instructions). As an alternative the parts can also be boiled (except the device).

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Qd 148-4_CE0124

ATMOS General Standard Terms and Conditions

1. General:

Our General Standard Terms and Conditions apply exclusively. Client's terms and conditions which are contrary to or deviate from our General Standard Terms and Conditions are not recognised unless their validity is explicitly confirmed in writing. Our General Standard Terms and Conditions also apply even if we deliver to clients without reservation, in the knowledge of the client's contrary terms and conditions. Our General Standard Terms and Conditions also apply to all future business with that client.

2. Proposal - Order Confirmation

Our proposals are subject to change without notice unless otherwise stated in our order confirmation. Each order is only accepted by us following our written order confirmation.

3. Orders

Every order requires an exact description of all of our product's details. We assume no liability for errors and damage caused by inaccurate or incomplete ordering details.

4. Prices

Unless otherwise stated in the order confirmation, our prices in the order confirmation are ex factory prices and exclude packaging and value added tax. Packaging is charged separately at cost price in the invoice. Value added tax is charged separately in the invoice according to the legal rate on the invoice date. We reserve the right to change prices appropriately should price reductions or increases, especially due to wage settlements, changes in the price of materials or currency fluctuations, be incurred. Proof of such changes will be provided for the client on request.

5. Payment Conditions - Balancing

Unless otherwise stated in the order confirmation, our invoices are payable with a 3% discount within 10 days (except for repair and assembly services) or within 21 days from the invoice date net cash; money receipts is decisive for complying with this term. We are entitled to charge interest after the due date at a rate 2% above the relevant basic interest rate of the German Federal Bank. Should the client have payment arrears, we are entitled to charge interest on arrears at a rate 5% above the relevant basic interest rate of the German Federal Bank. Should we be able to prove higher damages due to arrears, we are also entitled to claim these. The client only has the right to balance invoices against its own claims should such claims be confirmed in a court of law or recognised by us. The client does not have the right of retention due to disputed counterclaims.

6. Delivery Periods

Fulfillment of our delivery duties requires the punctual and proper fulfillment of the client's duties. The right to defense on the grounds of an unfulfilled contract is reserved. Should the client default in accepting the goods delivery or breach other cooperation duties, we are entitled either to withdraw from the contract or claim compensation for any increased costs incurred up to that time without setting a further deadline. The right to make further claims is reserved. Furthermore, in such cases, the risk of coincidental destruction or a coincidental deterioration in the quality of the delivered goods is transferred to the client in the case of default in accepting such goods or payment arrears. Acts of God or stoppages (due to insufficient supplies of material, industrial disputes etc.) entitle us either to demand an appropriate extension of delivery periods or to partly or entirely dissolve the delivery contract. This does not give the client the right to claim damages. We have fulfilled delivery periods if the delivery goods have left our factory or the client has been informed of the goods' readiness for delivery within such delivery periods. Delivery periods stipulated by the client are not recognised by us unless they form part of our order confirmation. We adhere to legal terms and conditions in cases where, as a result of an undue delay in the delivery for which we are liable, the client is entitled to claim that his interests in a continued fulfillment of the contract have ceased. We also adhere to legal terms and conditions should a delay in delivery be caused by deliberate or grossly negligent action by us or our representatives for which we are responsible. We are also responsible for such actions by our representatives or agents. Should the delivery delay not be caused by our deliberate infringement of contractual duties for which we are responsible, our liability is limited to damage which is regarded as typical for that case. We are liable according to the legal terms and conditions if and in so far as the delivery delay for which we are responsible is caused by an infringement of a substantial contractual duty. In such cases, our liability is also limited to damage which is regarded as typical for that case. Should the delivery delay be caused by a culpable infringement of non-substantial contractual duties, our client is also entitled

to claim a one-off damage compensation worth 3 percentage points of the delivery value of the goods for each week's delay, up to a maximum which is no higher than 15 percentage points of the delivery value of the goods

7. Delivery - Familiarisation

In the case of the delivery of devices for the medico-technical industry which require assembly and/or familiarisation for the final customer using specialist trade personnel (such as Ear, Nose and Throat Apparatus and Suction Units), we reserve the right to deliver the goods exclusively to the relevant specialist traders. Should the trader not carry out assembly and/or familiarisation for the final customer, this is carried out by us. In such cases, we reserve the right to charge the client for the additionally created costs. Our specialist traders operate a recording system so that, if necessary, our products can be traced to the final customer. The specialist trader undertakes to immediately report to us all events and risks which must be reported in connection with our products.

8. Passage of Risk - Packaging

Unless otherwise stated in our order confirmation, delivery is agreed ex factory. The risk of the goods' damage or loss is therefore transferred to the client as soon as the goods leave the factory or the client is in default of acceptance of the goods. This also applies to cases where we confirm prepaid carriage. Transport packaging and all other packaging according to the packaging regulations is not returnable. Our client is responsible for disposing of the packaging at its own cost. Our deliveries are insured by us at the client's expense unless explicitly otherwise agreed. No insurance is arranged in the case of goods which are collected by our clients. In the case of transport damage, claims are only handled if the client receives confirmation of any damage, reduced weight or loss by the shipping company before accepting the delivery.

9. Warranty

The client is responsible for examining the delivered goods immediately after receiving them to determine any eventual deficiencies or delivery errors, and to report these immediately. Should the client fulfil this examining and reporting responsibility, and should payment conditions be fulfilled, we shall be liable to the client within the scope of legal regulations. Our period of warranty shall in all cases be two years. Our client can make use of the warranty as follows, so long as he can provide first buyer proof (in the form of an invoice or delivery note) and provided that the product still has the original, unchanged serial number:

- We choose whether to fulfil our guarantee by providing repair services free of charge - either on the client's premises or in our factory - or replacing the product. We can also provide these guarantee services through an authorised company;
- Should a product be returned to us, the client agrees to send the product in its original or similar packaging, offering the same protection as the original packaging, to our address or any address notified by us.
- Our guarantee ceases to apply if changes of any kind have been made to our product, unless such changes have been made by us or a company authorised by us, or have been previously agreed upon in writing by us. Our guarantee also ceases to apply if third parties have carried out repairs to our products or replaced parts thereof. This applies regardless of the fact whether these measures individually or collectively led to a deficiency of the product;
- We accept no responsibility for damage defects caused by
 - operational wear and tear;
 - incorrect installation or incorrect or insufficient maintenance;
 - incorrect operation of the product (in contradiction to the handbook delivered with the product);
 - improper use or operating faults;
 - inappropriate or negligent handling and care, especially with respect to dirt, lime, suction of fluids, inappropriate cleaning and sterilisation;
 - using accessories and/or replacement parts which are not explicitly approved;
 - incorrect assembly and/or initial operation by the client or third parties;
 - the client's negligence in handling the product;
 - unacceptable operating conditions, such as humidity, temperatures, the power supply, vibrations.
 - accidents, acts of God, especially lightning, water, fire, public unrest and insufficient ventilation. We are not liable for damage to other objects apart from our product itself, except in the case of any deliberate or grossly negligent actions by us or our representatives or agents. Should no deliberate breach of contract be claimed, our liability is limited to damage which is regarded as typical for that case. This also applies in the case of our culpable infringement of substantial contractual duties. The indispensable conditions of German Liability Law remain unaffected thereby.

- For second-hand equipment, the period of warranty shall be reduced to a period of twelve months.

10. Reservation of Ownership

We retain ownership of our goods until the receipt of all payments arising from the business relationship, including all demands arising from installation orders, subsequent orders, repairs, accessory deliveries and replacement orders. Should we have agreed upon payment on the basis of cheque and bill transactions, the ownership reservation applies until the cheque received by us has been paid in, and does not expire through our credit upon receiving the client's cheque. In the case of a breach of contract by the client, especially payment arrears, we are entitled to repossess our goods. Repossession of our goods represents a withdrawal from the contract, unless explicitly declared in writing by us. We have the right to utilise the product after its repossession, whilst the income form such use is balanced against the client's arrears, after deducting appropriate utilisation costs. The client is responsible for handling the goods with care. Should maintenance and inspection work be necessary, the client must carry these out punctually at his own cost. Our client is entitled to sell the goods he has bought from us in a proper sale transaction. However, he must immediately assign all outstanding claims to the value of the final invoice sum (including value added tax) of our claims to his customers or third parties. The client is entitled to collect this claim even after such assignment. Our right to collect the claim ourselves remains unaffected thereby. We undertake to release the securities to which we are entitled if requested to do so by the client should the realisable value of the our securities be more than 10 percentage points higher than the outstanding claims. We reserve the right to choose the securities to be released.

11. Plans and Illustrations

We retain ownership of and copyrights to all plans, illustrations, calculations and other documents which are attached to our proposals. The client must receive explicit written permission before passing these on to third parties. Imitating our legally patented products is forbidden and will be prosecuted.

12. Jurisdiction and Place of Performance

Our central office is the place of performance for all disputes in connection with these General Standard Terms and Conditions and the contracts closed with clients under them. This jurisdiction excludes other jurisdiction relating to persons or subject-matter. Furthermore, our client is not entitled to bring charges against us in another court should he file counter-charges, carry out counterbalancing or declare retention. We, however, are entitled to bring charges against our client at their general place of jurisdiction or at another relevant court recognised by German or foreign law. Unless otherwise stated in the order confirmation, our central office is the place of performance.