**This addendum updates information in Trilogy100, Trilogy200, Trilogy202, and Trilogy O2 (all manuals).**

### Mouthpiece Ventilation (MPV)

MPV is a new feature in AC and PC modes for patients who are not invasively ventilated. The therapy, triggering, alarms, and user interface will provide on-demand ventilation without the need for an exhalation device. In addition, the Philips Respironics MPV Disposable Circuit further simplifies the user experience and ensures proper operation.

A new feature within MPV is the Circuit Disconnect (MPV) alarm. This alarm will notify the patient or caregiver that a spontaneous breath has not been initiated within the selected period. This may be useful to alert a caregiver when a patient is unable to connect to the mouthpiece or if the patient has fallen asleep.

### Setting up the MPV feature

#### Breathing Circuit Types

The available Philips Respironics MPV Disposable Circuit (or an equivalent circuit) is required for the proper operation of MPV. This circuit is supplied without an exhalation device and is not intended to support exhaling into the circuit.

For those patients who prefer to exhale into the breathing circuit, the passive and active circuits currently offered with the Trilogy ventilator may be used with MPV. Refer to chapter 4 of the manual (see section titled "Connect the Breathing Circuit to the Ventilator").

#### Connect the MPV Breathing Circuit to the Ventilator

When using the available Philips Respironics MPV Disposable Circuit, always follow the instructions included with the circuit.

#### Set up the Ventilator to Use MPV

1. Go to the Set-up screen when therapy is OFF. MPV cannot be enabled when the therapy is being delivered.
2. Select the Circuit Type: If you have a prescription that uses an active circuit, set the Circuit Type to match the prescription you normally use. Otherwise, set the Circuit Type to Passive.
3. Set the Mode to AC or PC. MPV is only available in these two modes.
4. Set MPV to ON.

5. Adjust the therapy parameters to the desired settings:
  - a. Tidal Volume (AC mode only): Range of 200 – 2000 ml
  - b. Breath Rate: Range of 0 – 60 BPM (Default: 0 to support the on-demand nature of MPV)
  - c. Inspiratory Time: Range of 0,3 – 5,0 seconds
  - d. Flow Pattern (AC mode only): Ramp or Square
  - e. IPAP (PC mode only): Range of 4 – 50 cmH<sub>2</sub>O
  - f. Rise Time (PC mode only): Range of 1 – 6
6. Adjust the alarm parameters to the desired settings:
  - a. Low Inspiratory Pressure alarm (AC mode only): Range of 1 – 40 cmH<sub>2</sub>O (Default: 1 cmH<sub>2</sub>O)
  - b. High Inspiratory Pressure alarm (AC mode only): Range of 10 – 80 cmH<sub>2</sub>O
  - c. Apnea alarm: Range of 10 – 60 seconds in 5 second increments (Default: 60 seconds). This alarm can be disabled by setting to OFF.
  - d. Circuit Disconnect (MPV) alarm: Available settings are: OFF, 1 minute, 2 minutes, 5 minutes, 10 minutes, 15 minutes (Default: 5 minutes).

## MPV Troubleshooting

**Question:** Why am I getting a Check Circuit alarm?

**Answer:** You will get a Check Circuit alarm when the breathing circuit is detached from the side of the ventilator. Re-attach the breathing circuit.

**Question:** Why am I getting a Low Inspiratory Pressure alarm when MPV is ON?

**Answer:** You will get a Low Inspiratory Pressure alarm when there is a large leak in the circuit. Check the circuit for loose connections, cracks, or holes.

**Question:** Why am I getting a Circuit Disconnect (MPV) alarm when MPV is ON?

**Answer:** You will get a Circuit Disconnect (MPV) alarm when you have not connected to the mouthpiece to initiate a breath within the time specified in the Circuit Disconnect (MPV) alarm setting. Check the setting and adjust as needed. The Circuit Disconnect (MPV) alarm can be disabled by the prescribing physician.

**Question:** Why am I getting extra air or an extra breath after I disconnect from the mouthpiece?

**Answer:** The Inspiratory Time setting may be too long or too short. Adjust the Inspiratory Time setting to match the time you are connected to the mouthpiece when you inhale. If the Inspiratory Time is too short, the breath will terminate early and the ventilator may deliver an extra breath after you pull away from the mouthpiece. If the Inspiratory Time is too long, the breath will continue after you have pulled away from the mouthpiece.

**Question:** Why am I getting an Apnea alarm when MPV is ON?

**Answer:** The Apnea alarm operates the same as it does in other modes, and occurs when you have not triggered a breath within the time specified in the Apnea alarm setting. Check the setting and adjust as needed. The Apnea alarm can be disabled by the prescribing physician.

## The Display and the MPV Feature

When MPV is ON, some of the patient parameters are not displayed. Only parameters that are applicable for the use of the ventilator with MPV are displayed:

- Pressure
- RR (Respiratory Rate)
- Vti (Inhaled Tidal Volume)
- Peak Flow
- PIP (Peak Inspiratory Pressure)

When MPV is ON, “MPV” will be displayed in the upper left corner of the monitoring screen.

*Note: When MPV is ON, the above patient parameters are measured at the outlet of the device.*

## AVAPS-AE Therapy

AVAPS-AE is a Bi-level therapy mode that automatically adjusts Expiratory Positive Airway Pressure (EPAP), Pressure Support, and the backup breath rate. AVAPS-AE monitors the resistance in the patient’s upper airway and adjusts EPAP automatically to maintain a patent airway. AVAPS-AE mode also monitors delivered tidal volumes and automatically adjusts pressure support to maintain the designated target tidal volume. AVAPS-AE also has the ability to automatically set and maintain a backup breath rate based on the patient’s own spontaneous breathing rate. AVAPS-AE is intended for non-invasive applications only.

AVAPS-AE is enabled in the Settings and Alarms Screen. When in AVAPS-AE mode, the AVAPS feature is always enabled. Similar to AVAPS, AVAPS-AE supports tidal volumes from 50 to 2000 ml (in 5 ml increments).

The following new parameters should be set to ensure appropriate therapy delivery and maximize patient comfort:

- **Maximum Pressure:** When in AVAPS-AE mode, the Maximum Pressure delivered can be set from 6 to 50 cmH<sub>2</sub>O in increments of 1. Settings lower than the maximum setting of 50 cmH<sub>2</sub>O may limit other pressure settings.
- **Pressure Support Max:** When in AVAPS-AE mode, Pressure Support Max can be set from 2 to 40 cmH<sub>2</sub>O in increments of 1. This setting may be limited by the Maximum Pressure setting and EPAP Min Pressure setting.
- **Pressure Support Min:** When in AVAPS-AE mode the minimum pressure support may be set from 2 to 40 cmH<sub>2</sub>O in increments of 1. It may not exceed the Pressure Support Max setting. This setting may be limited by the Maximum Pressure or EPAP Max Pressure setting.
- **EPAP Max Pressure:** When in AVAPS-AE mode, the EPAP Max Pressure can be set from 4 to 25 cmH<sub>2</sub>O in increments of 1. This setting cannot be lower than the EPAP Min Pressure and may be limited by the Maximum Pressure setting and Pressure Support Min setting.
- **EPAP Min Pressure:** When in AVAPS-AE mode, the EPAP Min Pressure can be set from 4 to 25 cmH<sub>2</sub>O in increments of 1. This setting cannot exceed EPAP Max Pressure and may be limited by the Maximum Pressure setting and Pressure Support Min setting.
- **Breath Rate:** When in AVAPS-AE mode, the breath rate can be set from 0 (off) to 60 BPM. The breath rate can also be set to “Auto”.

*Note: When “Auto” Breath Rate is selected, the Inspiratory Time will be set by the device and is not user settable.*

## AVAPS-AE Troubleshooting

**Question:** Why can't I increase my Pressure Support Max to the full range of 40 cmH<sub>2</sub>O?

**Answer:** Some settings, such as Pressure Support Max, can be limited by other pressure settings. For example, reducing the Maximum Pressure setting will impose the same reduction on the Pressure Support Max capability of the device. Increasing the EPAP Min Pressure will also reduce the Pressure Support Max that can be applied. You may have to modify one setting in order to extend the range of another setting.

**Question:** What happens if the device needs to increase EPAP but the Maximum Pressure setting has already been reached?

**Answer:** The EPAP will increase as high as the EPAP Max Pressure setting. The amount of Pressure Support applied may be reduced in order to reach the EPAP Max setting. Airway patency and hence the EPAP pressure will take precedence over Pressure Support in these situations.

**Question:** How does the Ramp function work when using AVAPS-AE?

**Answer:** Ramp: Activation of the Ramp Feature in AVAPS-AE will reset the EPAP setting to the EPAP Min Pressure and restart the Auto EPAP algorithm. Additionally, the Pressure Support will be reduced to the Pressure Support Min level and ramp as high as the Pressure Support Max, as required by the AVAPS algorithm, to allow for a comfortable transition to sleep.

*Note: If using AVAPS, the pressure support capability will ramp from the IPAP Min (or Pressure Support Min) to the IPAP Max (or Pressure Support Max).*

### The AVAPS feature has been changed to allow for the adjustment of the rate.

The AVAPS Rate setting allows you to adjust the maximum rate at which the pressure support automatically changes to achieve the target tidal volume. The actual rate may be less than this maximum setting depending on how far the current estimated tidal volume is from the target tidal volume. A higher rate allows the AVAPS algorithm to change pressure support faster to meet the target tidal volume. It can be set from 1 cmH<sub>2</sub>O per minute to 5 cmH<sub>2</sub>O per minute in increments of 1 cmH<sub>2</sub>O per minute.

### The following warning and associated notes are removed in all Trilogy manuals.

C-Flex, Bi-Flex, and AVAPS are intended for use by adult patients.

### The following information replaces SIMV and PC-SIMV therapy modes in Chapter 3 under Therapy Modes.

SIMV/PC-SIMV mode provides a mixture of mandatory and spontaneous breaths. SIMV/PC-SIMV mode will guarantee one mandatory breath in each SIMV/PC-SIMV cycle window. This mode can be used with either volume or pressure targeted ventilation. Spontaneous breaths can be delivered with pressure support if selected.

The SIMV/PC-SIMV breath rate determines the length of the SIMV/PC-SIMV cycle window. The first phase of the SIMV/PC-SIMV cycle window is reserved for synchronizing a mandatory breath with patient effort. If the patient triggers a breath during this phase of the cycle window, the ventilator will deliver a synchronized mandatory breath. Any breaths triggered after a mandatory breath in any SIMV/PC-SIMV cycle window will be spontaneous breaths. The ventilator will automatically deliver a mandatory breath if a patient fails to trigger a breath during the mandatory phase of the SIMV/PC-SIMV cycle window. This process will be repeated at the start of every SIMV/PC-SIMV cycle window.

## The Trilogy display has been revised for ease in distinguishing alarm parameters from therapy settings.

The alarm parameters have been renamed and will be displayed in the Settings and Alarms menu. “Alarm” has been added to each parameter name. For example, “Circuit Disconnect” is now displayed as “Circuit Disconnect Alarm”. The processing and description of the alarms has not changed.

## A new setting has been added to increase the sensitivity of Auto-Trak.

Sensitive Auto-Trak is an enhancement to the Auto-Trak algorithm that improves patient and machine synchrony for patients with minimal respiratory effort. Sensitive Auto-Trak refines the baseline trigger sensitivity threshold. The device can be set to trigger breaths based on automatic flow thresholds or specific flow settings. Select Auto-Trak [Sensitive], Auto-Trak, or Flow Trigger as the Trigger Type. When Auto-Trak or Auto-Trak [Sensitive] is selected, the trigger initiates based on automatic flow thresholds. When Trigger Type is set to Flow Trigger, Flow Trigger Sensitivity, Leak Compensation, and Flow Cycle Sensitivity become active, and the trigger initiates based on the Flow Trigger Sensitivity setting.

## The following information replaces Low Vte Alarm in Chapter 3 under Ventilator Alarms.

This is a high priority alarm. It occurs when the estimated exhaled tidal volume is lower than the Low Vte alarm setting for three consecutive breaths. The device continues to operate. The alarm will automatically terminate when a breath occurs in which the exhaled tidal volume exceeds the Low Vte alarm setting.



When AVAPS is on, this alarm will occur when the calculated tidal volume is less than or equal to 90% of the target tidal volume setting for at least 1 minute. The alarm will automatically terminate when the calculated tidal volume exceeds 90% of the target tidal volume for at least 30 seconds.

**Note:** The Pressure Support Max limit must be reached to enable this alarm.

## Clinicians may view pressure and flow in waveform format by using the Waveform screen.

The waveforms are only available in Full Menu Access mode. In order to access the Waveform screen, select the “Waveform View” item in the main menu. In addition to pressure and flow waveforms, numerical values are also displayed for respiratory rate, tidal volume, and leak.

The Waveform screen can be adjusted in order to improve viewing. The Right (Graph) button changes the functionality of the Up and Down buttons each time it is selected among the following options:


- **Pause/Play:** When the Waveform screen is initially displayed, the Up and Down keys will perform the Pause (  ) function. When either the Up or Down button is selected, the graphs will stop updating and the functionality of the Up/Down buttons will change to Play (  ). When the Up or Down buttons are selected for Play, or there is no user interaction within one minute of activating the Pause function, the Up/Down buttons change back to Pause, the graph data will be erased, and the graphs will begin to update from the left side of the graphs.
- **Time Scale:** Selecting the Up and Down keys while the Time Scale (  ) function is active adjusts the scaling of the X-axis for the graphs. Pressing either button will cycle through the available time scales and adjust the currently displayed data for the new time scale.
- **Data Scale:** Selecting the Up and Down keys while the Data Scale (  ) function is active adjusts the Y-axis of the graphs. Pressing either button will change the scaling of all graphs to best fit the current data.

Selecting the Left (Exit) button will return you to the Monitor screen. If an alarm is generated while the Waveform screen is displayed, the ventilator will automatically exit the Waveform screen and display the Alarms and Messages Screen where the details of the alarm will be displayed. The Monitor screen will also automatically be displayed if there is no user interaction for 10 minutes when the Waveform screen is displayed.


### Trilogy allows the addition of nebulizer treatments without triggering unnecessary alarms.

During Nebulizer Treatment with a Passive circuit, the Low Circuit Leak alarm sensitivity is adjusted for a 20-minute period in order to reduce nuisance alarms. The Low Circuit Leak alarm will still be sensitive enough to alert the user of potential hazards, such as a blocked exhalation port, but the number of nuisance alarms should be significantly reduced.

The Nebulizer Treatment feature can be enabled from the Settings and Alarms menu when using a Passive Circuit. When enabled, a Nebulizer on-screen button appears on the Monitor screen. When pressed, the Nebulizer Treatment period starts. Once the 20-minute Nebulizer Treatment period has expired, the Low Circuit Leak alarm sensitivity will return to its normal setting. Once the Nebulizer Treatment period is started, the Cancel Nebulizer button and the Nebulizer icon will be displayed. Pressing the Cancel Nebulizer button will stop the 20-minute window, remove the Nebulizer icon, change the button back to the Nebulizer button, and return the Low Circuit Leak alarm sensitivity to its normal setting.

Nebulizer Enabled Screen

Nebulizer Button Screen

### Trilogy now allows for pressure support up to 40 cmH<sub>2</sub>O.

IPAP, IPAP Max, and IPAP Min can now be set as high as 40 cmH<sub>2</sub>O above EPAP. In SIMV/PC-SIMV mode, the Pressure Support and/or Inspiratory Pressure can be set as high as 40 cmH<sub>2</sub>O above PEEP. The maximum pressure delivery of the Trilogy remains at 50 cmH<sub>2</sub>O.

### The following information replaces a portion of the instructions in Chapter 7 under Cleaning the Ventilator.

Use a cloth with one of the following cleaning agents to clean the exterior of the device.

- Mild detergent
- 70% Isopropyl Alcohol
- DisCide Towelettes
- 10% Chlorine bleach solution

### The following information has been added to the Standards Compliance listing in Chapter 11.

- IEC 60601-1-8: General Requirements for safety – Collateral Standard: Alarm systems in medical electrical equipment and medical electrical systems
- IEC 60601-1-11: General Requirements for Safety and Essential Performance – Collateral Standard: Medical equipment and medical electrical systems used in the home healthcare environment

The degree of ingress protection for the Trilogy device has changed. The following information updates the Electrical specifications in Chapter 11.

Degree of Ingress Protection: Device Drip Proof

IP22

The following warning has been added to the Electro Magnetic Compatibility (EMC) information.

The use of accessories, transducers and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

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Respironics Inc.