PERCUSSIONAIRE®

Travel Air® User Manual

TRUE-IPV® Therapy Device

EN

AR

DE

ES

FR

HE

IT

JA

NL

NO

RU

Dear Customer,

Congratulations on the purchase of your Travel Air®! Your new TRUE-IPV® therapy device is portable, self-contained, and easy to use. It allows you the freedom to go where you want to go while complying with the treatment plan prescribed by your doctor!

To make sure you receive maximum benefit from TRUE-IPV® therapy and to ensure your safety, please read the enclosed Travel Air® User Manual. It contains important SAFETY and TECHNICAL DATA that should be kept handy for easy reference. Your safety and satisfaction are important to us!

For assistance setting up or using the Travel Air®, or to report unexpected operation or events, please contact Percussionaire® or your local distributor. With proper care, the Travel Air® will reward you with long, trouble-free service life.

Thank you for placing your trust with us!

Sincerely,

Percussionaire® Corporation

© 2020 Percussionaire® Corporation

ALL RIGHTS RESERVED

1st Edition

First Printing September 2020

Percussionaire® is a registered trademark of Percussionaire Corporation.

This work is the sole property of Percussionaire® Corporation. The information in this manual is confidential and may not be disclosed to third parties without the prior written consent of Percussionaire®. No part of this document may be copied, reproduced, transmitted, or stored in any electronic information system without the prior written consent of Percussionaire® Corporation.

Authorized European Union Representative:

MDSS GmbH, Schiffgraben 41, 30175 Hannover, Germany

The devices and products contained in this manual may be covered by one or more patents.

This manual was originally released and supplied in English.

For a list of available translations, contact customerservice@percussionaire.com.

All ventilators should be operated and serviced only by trained professionals. Percussionaire® Corporation's sole responsibility with respect to its ventilators, accessories, components, and software, and their use, is as stated in the warranty provided in the manuals. The information set forth herein is believed to be accurate; it is not a substitute for the exercise of professional judgment.

Table of Contents

Chapter 1: Introduction	1
Intrapulmonary Percussive Ventilation (IPV®)	1
Chapter 2: Intended Use	2
Indications for Use	2
Patient Populations	2
Absolute Contraindications	2
Relative Contraindications.....	2
Possible Adverse Reactions.....	2
Physiological Benefits of TRUE-IPV®	2
Clinical Limitations/Restrictions	3
Document Symbols	3
Chapter 3: System Description	4
Front Panel	4
Rear Panel	4
Function Overview	5
Digital Multimeter (PDM)	6
Power-On Self-Test (POST) Mode	6
Wake Mode	6
Active Mode	7
Report Mode.....	7
Sleep Mode	8
Fault Mode	8
Fault Logging.....	8
Fault Detection	9
Phasitron® P5-TA Breathing Circuit	9
Phasitron® P5-TA Diagram	10
Configurations	10
Chapter 4: Travel Air® Setup	11
Install the Phasitron® P5-TA Holder Bracket	11
Attach the Handle	11
Connect the Breathing Circuit Tubing	12
Connect to the Power Supply	12
Digital Multimeter (PDM) Setup	13
Removal of PDM from device	13
Connect to the Phasitron® P5-TA	14
Removing Nebulizer Bowl	14
Adding Saline, Sterile Water, or Medication	15
Reattaching Nebulizer Bowl.....	15
⚠ CAUTION:.....	15
Timestrip®	16
Attach Timestrip® to Phasitron® P5-TA Bag.....	16
Chapter 5: Travel Air® Pre-Use Check	17

Chapter 6: TRUE-IPV® Treatment	19
Before Starting Treatment	19
Amplitude	19
Percussion Frequency.....	19
Pulse Ratio	19
Getting Ready	19
Travel Air® Controller	20
Starting Treatment	20
Chapter 7: Travel Air® Indicators	22
High-Pressure Indicator Alarm.....	22
High-Temperature Shutdown	22
Chapter 8: Cleaning and Disinfection	23
Homecare Cleaning Schedule.....	23
Travel Air® Controller	23
Phasitron® P5-TA Breathing Circuit Kit.....	24
Disassembly of the Phasitron® P5-TA	24
Cleaning Procedure	25
Disinfection Solution Procedure	26
Dry Inside of Tubing.....	27
Phasitron® P5-TA Lubrication After Cleaning and Disinfection.....	27
Reassemble the Phasitron® P5-TA	28
Chapter 9: Maintenance	30
Home Maintenance Schedule.....	30
Lubrication	30
Drain Water Trap.....	31
Filter Replacement	31
Water Trap Gasket	32
Patient Air Filter Replacement.....	32
Cooling Air Filter	33
Fuses	34
Switch Fuse	34
Thermal Fuse	34
PDM (Multimeter) Batteries.....	35
Changing PDM Batteries	35
Chapter 10: Accessories	36
Chapter 11: Troubleshooting	36
Chapter 12: Service and Repair	37
Disposal of Equipment.....	37
Returns.....	37
Chapter 13: Technical Specifications	38
Travel Air® Controller	38
Phasitron® P5-TA Specifications	39
Multimeter (PDM) Specifications.	39
EMC Warnings	40
Guidance and Manufacturer’s Declaration – Electromagnetic Emissions	40
Recommended Separation Distances	43
Chapter 14: Limited Warranty	44
Chapter 15: Label	44

Chapter 1: Introduction


This chapter provides an overview of the Travel-Air® device and TRUE-IPV® therapy.

Intrapulmonary Percussive Ventilation (IPV®)

The Travel Air® Impulsator® is a portable, self-contained, easy to use device that allows patients to benefit from TRUE-IPV® therapy at home. Designed specifically for non-continuous home care use, the Travel Air® is a pressure-limited, flow-interrupted, and time-cycled ventilator that delivers TRUE-IPV® therapy, a modality of mechanical ventilation.

The Travel Air® device delivers high-frequency percussive pulses between 60-500 pulses per minute. The high-frequency pulses ramify throughout the airways and alveolar ducts, augmenting diffusive ventilation in the gas exchange regions of the lungs, allowing improved FRC, CO₂ removal, airway clearance, and lung recruitment.

During operation, the Travel Air® system supplies a continuous, dense aerosol mist of saline or sterile water, which is delivered into the lungs during therapeutic percussion, serving to reduce the adhesive and cohesive forces of retained airway secretions.

The Travel Air®, using a Phasitron® P5-TA, provides intrapulmonary percussive ventilation either invasively through an artificial airway or non-invasively, by mouthpiece or mask. The patient can turn the device On and Off and adjust the amplitude and frequency.

Chapter 2: Intended Use

Indications for Use

The Travel Air® is indicated for the mobilization of secretions, lung expansion therapy, and the treatment and prevention of pulmonary atelectasis.

Patient Population

The Travel Air® Impulsator® is for use on patient populations older than 2 years through adult.

Absolute Contraindications

Not for use as a critical care, life-support device.

- | | |
|----------------------------------|-----------------------------------|
| • Untreated tension pneumothorax | • Untrained or unskilled operator |
|----------------------------------|-----------------------------------|

Relative Contraindications

• History of pneumothorax	• Lack of patient cooperation
• Recent pneumonectomy	• Vomiting
• Pulmonary hemorrhage	• Pulmonary air leak (without functioning chest tube)
• Myocardial infarction	

Possible Adverse Reactions

• Decreased cardiac output	• Increased intracranial pressure
• Pneumothorax	• Increased air trapping
• Hyper-oxygenation	• Pulmonary air leak
• Pulmonary hemorrhage	• Hyperventilation
• Gastric distension	

Physiological Benefits of TRUE-IPV®

• Recruitment of atelectatic lung	• Mechanical bronchodilation
• Improved FRC	• May improve breathing pattern
• Decreased work of breathing	• Increased secretion mobilization

Clinical Limitations/Restrictions

Use of the Travel Air® is limited to individuals who have received proper training.

 **WARNING:** Because pulmonary alveoli cannot be ventilated when their transmitting airways are obstructed, suctioning should be performed as necessary.

NOTE: A **WARNING** icon indicates a risk of injury to patient or operator.
A **CAUTION** icon indicates a risk of equipment damage.

Document Symbols

 WARNING	 Type BF Applied Part
 CAUTION	 Single Patient Use
 Read the manual before use	 Prescription Only
 CE marking	 Catalog Number
 Manufacturer	 Lot Number
 Manufacture Date	 European Representative
 Non-Sterile	 Not Made with Natural Rubber Latex
 Does Not Contain the Phthalate Plasticizers DEHP, DIBP, DBP, or BBP	 Disposal

EN


Chapter 3: System Description

Front Panel

1. Multimeter display	2. Breathing circuit connections
3. Power indicating LED lights	4. Percussion/Frequency knob
5. Amplitude knob	

Rear Panel

6. Power switch	7. Power switch fuse
8. Thermal fuse	9. Release latch for cooling air filter
10. Power cord socket	11. Handle mount
12. Patient air inlet filter behind water trap bowl	13. Cooling air outlet
14. Water trap	15. Cooling filter air inlet

Function Overview

Digital Multimeter (PDM)

The Percussionaire® Digital Multimeter (PDM) has six different operating modes: POST, Wake, Active, Report, Sleep, and Fault.

Percussion

The **PERCUSSION** control knob adjusts the frequency of percussions distributed to the patient every minute. The number of percussions ranges from approximately 60 to 500 pulses/minute (1 to 8.3 Hz).

Amplitude

The amplitude control knob adjusts the volume delivered to the patient.

NOTE: *Amplitude is affected by lung compliance and resistance. Amplitude is attenuated by resistance and compliance to an airway pressure.*

Breathing Circuit Tubing Sockets

The Phasitron® P5-TA breathing circuit tubing connectors are inserted in the like-colored tubing sockets for easy setup.

Power Switch

Located on the back of the device, the power switch turns the Travel Air® ventilator on and off.

EN

Digital Multimeter (PDM)

The Digital Multimeter (PDM) has six different operating modes: POST, Wake, Active, Report, Sleep and Fault

Power-On Self-Test (POST) Mode


When batteries are installed in a system, the Percussionaire® Digital Multimeter (PDM) software displays the software revision, battery voltage, total usage time, and serial number for 15 seconds. This Start-Up mode allows the software to perform additional tests on the hardware that are part of the **Power-On Self-Test (POST)**. If any errors are detected in POST, the PDM enters the Fault mode.

NOTE: Do not install PDM until the POST check is complete, and the screen is blank, indicating Sleep mode.

System Information Display

```
Percussionaire  
Digital Multimeter  
(C) 2014, RDI  
Bat: 3.05V  
Total Time: 23,075h 27  
Code Rev: 2.XX  
Serial #: 2140604-001
```

Rear view of PDM

Measuring port connection
Do NOT touch!


Reset button
(only used
for firmware
upgrades)

Wake Mode

To wake up the PDM, ensure the ventilator pressure is greater than 2.5 cmH₂O/hPa at the Phasitron® P5-TA patient delivery port for more than 1 second.

The PDM remains on for the first 15 seconds, showing the Bar-Graph timer.

If usage is stopped within 12 seconds, the PDM enters Report mode. After 15 seconds, the current session continues counting from 16 seconds, which turns into Active mode.

The PDM display screen in Wake mode


NOTE: Display numbers are for reference only.

Model: Home IPV

Device: Travel Air® Impulsator®

Display Metrics: Pulse Frequency Rate, Pulse Amplitude Pressure

Active Mode

At 16 seconds the PDM enters Active mode. The timer bar will change to a numeric display, showing the current usage Session Timer. Above the timer is the pulse amplitude display. This is calculated from the pressure measurements at the moment of instantaneous peak and trough amplitude averaged over 5 seconds. The Active mode display also shows the currently measured percussion rate/pulse frequency.

The PDM displays the usage Session Timer in minutes and seconds. The Session Timer is the total time of the current usage. The Session Timer can display a maximum of 59 minutes and 59 seconds. If usage has been stopped for more than 5 minutes, the Session Timer will reset and start over.

NOTE: To display the most recent usage duration time, see Report mode.

Report Mode

A.		<p>The Session Timer and the Total Usage Timer (A) are displayed for 2 seconds, followed by the System Information page (B) for 2 seconds, alternating. The alternating page display continues for 5 minutes, or until usage resumes and the PDM enters Active mode.</p>
B.		<p>During the 5-minute period, a horizontal bar graph indicates the time by moving from left to right at a fixed rate. After 5 minutes of no usage, the system information page is no longer displayed, and the time display flashes (2 seconds on, 2 seconds off) (C) for an additional 25 minutes.</p>
C.		
		<p>When the Travel Air® is switched off, the measurements drop to zero after a few seconds. The PDM enters Sleep mode after 25 minutes.</p>


Sleep Mode

Blank Screen – PDM Sleep Mode

In **Sleep** mode, the LCD is off, but the microcontroller continues to sample and calculate the pressure at the measuring port 5 times a second. Over any 3-second period, if the pressure is greater than 2.5 cmH₂O/hPa at the Phasitron® P5-TA patient delivery port, for more than 1 second, the PDM enters the Wake mode.

Fault Mode

The PDM displays an error message on the LCD stating, “Contact Factory for Service” and stays in this mode until both batteries are removed. The displayed information includes the software revision, PDM serial number, the total usage time, and an error code for the exclusive use of the factory.

In all other modes, the software continuously monitors the hardware for errors and verifies that each data sample has a valid value. If an error is detected, the software logs the error and reboots the processor. Rebooting causes it to recover from a transient error. The processor then returns to the same mode it was in before the reboot. If more than one error is detected in any 10-second period, it is considered a fatal error, and the software enters **Fault** mode.

NOTE: Pressure faults are triggered by a continuous pressure of more than 150 cmH₂O for more than 5 seconds during Wake and Active modes.

NOTE: If the **System Failure** screen is displayed, remove batteries for 30 seconds. Replace batteries (note that positive terminals face the same direction) and wait 30 seconds until the screen turns off. If POST check runs correctly, PDM may be used. If the **System Failure** screen recurs, contact an authorized Percussionaire® service center.

Fault Logging

The software keeps track of several types of hardware and data faults. All faults are logged in the microcontroller’s memory and are retained even if the batteries are removed. If multiple faults happen within 10 seconds of each other, the PDM stops normal operation and enters Fault mode. In this mode, a subset of the collected fault information is displayed on the LCD. This data is intended for manufacturing and repair use only.

The user can exit the Fault mode by removing and replacing the batteries. This resumes normal operation of the PDM but does not erase the faults stored in memory or fix the problem that caused the fault.

Fault Detection

The PDM has both hardware and software fault detection. This is a dedicated hardware “watchdog” that runs on an independent clock source and can continue to operate even if the main microprocessor’s clock fails or the microcontroller pauses in any way. The independent fault detection is reset each time a valid pressure reading (free of hardware and software errors) is obtained.

In addition to the hardware fault detection, the software also implements a fault detection watchdog. This watchdog detects if a software task fails to complete within the specified time, logs an error, and resets the processor.

Phasitron® P5-TA Breathing Circuit

The patented Phasitron® P5-TA uses a unique sliding venturi mechanism to protect the lung from overpressure. By automatically adjusting to the resistance of the lung, the Phasitron® P5-TA precisely and safely delivers the optimal amount and pressure of air required by the alveolar space. When lung resistance is low, as in a compliant lung, all the pulsed air from the Travel Air® enters the mouth of the venturi. Each air pulse draws up to four times as much additional air into the venturi tube. This low-pressure entrained air automatically fills the available space in the lung. The Phasitron® P5-TA continuously and instantaneously adjusts to keep a gentle and safe air pressure, even in a compromised lung.

EN


NOTE: TRUE-IPV® therapy can only be achieved using the Phasitron® P5-TA.

⚠ WARNING: Read and understand the instructions completely before setting up and using your Phasitron® P5-TA breathing circuit kit.

⚠ WARNING: All persons performing home IPV® treatment must be trained in the use of the Travel Air® device.


EN-9

Phasitron® P5-TA Diagram

Configurations

Phasitron® P5-TA kit can be used with or without a mouthpiece or standard mask (as shown below). Connection sizes: 15 mm ID or 22 mm OD.

Chapter 4: Travel Air® Setup

Install the Phasitron® P5-TA Holder Bracket

Insert the bracket in the left or right slot, depending upon preference.


Attach the Handle

NOTE: *The handle is for transporting only. Remove the handle before use.*


Insert the handle into the mounting locations on the back of the Travel Air®. With the handle pressed firmly into the mounting locations, lift until a “click” is heard. To remove, retract the latch by pulling back on either side of the release at the location shown.

Connect the Breathing Circuit Tubing

Connect to the Power Supply

The power switch must be in the off position. Connect the power cord to the socket at the back of the Travel Air® and the electrical outlet.

NOTE: *The power cord can be used as the main disconnect for the Travel Air®.*

 **WARNING:** Make sure that the power cord is not damaged before connecting it to the Travel Air®.


Digital Multimeter (PDM) Setup

NOTE: To ensure correct atmospheric pressure calibration at startup, remove pull tabs or batteries, wait 30 seconds, and reinstall. Allow 15 seconds for power-on self-test. When the screen goes blank, install the multimeter into the device.


Remove PDM from Travel Air®

1. Press in on the PDM's bezel and twist counterclockwise approximately 20 degrees.
2. Gently pull on the PDM to remove it from the housing.
3. Install the PDM back into the housing and twist clockwise until the stop is felt.


NOTE: To access battery pull tabs, turn the PDM counterclockwise and remove from the device.


Side view of PDM


Connect to the Phasitron® P5-TA

1. Connect the yellow quick-connect fitting to nebulizer bowl.


2. Connect the red tubing onto conical connector at the rear of the Phasitron® P5-TA body.


3. Connect the clear quick-connect fitting to the cap at the rear of the Phasitron® P5-TA body.

Removing Nebulizer Bowl

1. The nebulizer bowl can be removed by rotating bowl towards the rear of the Phasitron® P5-TA.

2. Rotate gently past the stop and the nebulizer bowl will release.

Adding Saline, Sterile Water, or Medication


3. Twist clockwise to open nebulizer bowl. Add prescribed amount.


4. Reverse to close.

Reattaching Nebulizer Bowl

5. Attach past stop and rotate bowl towards front of the Phasitron® P5-TA.

6. Once bowl is in position, the Phasitron® P5-TA is ready for use.

⚠ CAUTION

⚠ **CAUTION:** Ensure the yellow nebulizer tubing is not bent. This may cause undue stress on connector.

⚠ **CAUTION:** Do not bend nebulizer bowl while holding the tubing. This may cause undue stress on the red line conical connector.

EN

Timestrip®

Bubble

Using the Timestrip®

It is very important to replace your Phasitron® P5-TA kit every 6 months for the most effective treatments. Timestrip® is a visual indicator that makes it easy for you to remember when to replace your Phasitron® P5-TA.

Activating Timestrip®

Firmly squeeze the bubble between your finger and thumb. You may need to squeeze it more than once to ensure all the liquid is emptied from the bubble.

Verifying Timestrip®

A red line will appear on the Timestrip® within 1 to 2 minutes, indicating activation. If you do not see the red line after 3 minutes, press the bubble again.

NOTE: *The Timestrip® is only supplied with the Phasitron® P5-TA home care kit.*

Attach Timestrip® to Phasitron® P5-TA Bag

Remove the adhesive tab on the back of the Timestrip® and attach Timestrip® to the front of your Phasitron® P5-TA storage bag.

When the red bar reaches 6 months, it is time to replace your Phasitron® P5-TA breathing circuit kit. Contact your home healthcare provider to order your new Phasitron® P5-TA kit.

Timestrip® does not monitor the condition of the Phasitron® P5-TA kit. Timestrip® only records the passage of time after it is activated. If activated at the time of Phasitron® P5-TA first use, Timestrip® serves as a reminder that the kit needs to be replaced.

The liquid inside Timestrip® is not toxic. If tampered with, the Timestrip® may leak and will not function properly.


⚠ WARNING: The Timestrip® may present a choking hazard to small children. Adult supervision is required.

⚠ WARNING: If the liquid inside the Timestrip® contacts skin, wash with soap and water. If the liquid contacts eyes, rinse immediately with water for 15 minutes. If irritation develops, seek medical attention.

Chapter 5: Travel Air® Pre-Use Check

Regularly perform the pre-use check before using the Travel Air® for treatment. If the Travel Air® has been stored in conditions outside of the operational conditions outlined in Chapter 13 of this manual, allow the unit to acclimate for 2 hours before using.

1. Check that the Travel Air® is clean on the outside and that the power cord is in good condition.
2. Attach the Phasitron® P5-TA and tubing correctly, connecting it to the Travel Air®. Place the Phasitron® P5-TA in the holder bracket.
3. Connect the power cord to the Travel Air® device.
4. Turn on the power and check that the LED lights on the front panel turn on.
5. Occlude the Phasitron® P5-TA patient port with the supplied cap or use finger.

6. Rotate the **PERCUSSION** (frequency) control knob full clockwise (right) to stop. Do not force past the natural stop.


7. Rotate the (left) Amplitude knob counterclockwise to between 90 and 99 cmH₂O/hPa without alarm activating.

8. Plug or pinch off clear line going to Travel Air and confirm alarm sounds.

9. Rotate the (right) **PERCUSSION** knob full counterclockwise until stop is felt. Rotate the (left) Amplitude knob full counterclockwise. With both knobs turned full counterclockwise, observe a frequency between 500 and 600 pulses/minute and a pressure between 30 and 90 cmH₂O/hPa.

10. Rotate the (left) Amplitude knob slowly, full clockwise. Do not force past natural stop. Observe a gentle decrease in flow until the device is off.

11. Rotate both knobs so that "+" symbols are at the 9:00 position and "-" symbols are at 3:00 position.

12. If all steps are performed and function as described, this completes the start-up test.

Chapter 6: TRUE-IPV® Treatment

A first treatment should always be low amplitude and high percussion (frequency), between 300 to 500 pulses per minute. The handheld Phasitron® P5-TA delivers TRUE-IPV® therapy as you breathe through the mouthpiece. You can adjust percussion and amplitude by using the simple controls on the Travel Air® device.

Before Starting Treatment

The patient should be trained in the proper use of home TRUE-IPV® therapy. Trained patients may adjust the amplitude and percussion settings under the guidance of their prescribing physician, practitioner, or clinician.

Amplitude

The (left) Amplitude control knob adjusts the amplitude of pulses to the Phasitron® P5-TA. The Phasitron® P5-TA functions as a safety clutch that entrains or diverts inspiratory air flow to the patient, depending on airway compliance and resistance.

Percussion Frequency

The (right) **PERCUSSION** knob controls the number of percussions distributed to you every minute. The number of percussions ranges from approximately 60 to 500 pulses/minute (1 to 8.3 Hz).

Pulse Ratio

The ratio of the duration of the pressure wave to the cycle is automatically adjusted by the **PERCUSSION** control knob.

Getting Ready

Patients should be able to close their lips tightly around the mouthpiece. A full mask can be used if the patient is unable to do this.


Phasitron® P5-TA
with mouthpiece

Phasitron® P5-TA
with face mask

Travel Air® Controller

1. Turn on the Travel Air® power switch located on the back panel. The LED lights will come on, and the compressor will start within 5 seconds.
2. Check that the Travel Air® is operating properly per the pre-use check procedure.
3. Rotate the **PERCUSSION** control knob counterclockwise to start the percussions.

NOTE: Frequency is not indicated on the multimeter display until treatment begins or the Phasitron® P5-TA is occluded at the mouthpiece. Verify that the tubing harness is correctly connected to the Travel Air® and the Phasitron® P5-TA.

Starting Treatment

Patients should sit in an upright position with the shoulders relaxed. An ordinary chair, with or without armrests, can be used. The treatment can be performed on patients in wheelchairs or lying down. Patients should breathe normally through the mouth (not through the nose).

While the following steps are descriptive of a typical TRUE-IPV® treatment, refer to the guidelines provided to you by your physician/clinician/respiratory therapist for your specific treatment.

1. Turn the (right) **PERCUSSION**/Frequency control knob full counterclockwise. This is where you will start.
2. Turn the (left) Amplitude control knob full clockwise to the off position, for nebulizer only.
3. Slowly rotate the Amplitude control knob counterclockwise until movement of the chest near the diaphragm is observed. Adjust according to your comfort level.

4. Turn the **PERCUSSION** (frequency) control knob clockwise to slow frequency and increase the amplitude to a comfortable range. Intrapulmonary Percussive Ventilation (IPV®) is generally more effective when it includes stronger pulsatile flow.

5. Continue treatment for 20 minutes, or as directed by the prescription.

NOTE: *IPV® treatment generates a considerable mobilization of secretions. Drainage techniques (such as controlled coughing) are particularly important for patients without spontaneous coughing.*

 **WARNING:** Avoid, if possible, treatment for at least one hour after eating or before going to bed.

6. When treatment is complete, turn off the power switch.

Chapter 7: Travel Air® Indicators

High-Pressure Indicator Alarm

Obstruction of the expiratory port by mucus or other blockages in the venturi can prevent return flow. This can increase the patient's expiratory pressure to 100 cmH₂O. If this happens, the patient will be protected by an audible indicator. In the event of a high-pressure indication – look for blockages or kinks in the circuit and tubing, beginning from the patient and working toward the Travel Air®. Remove any blockages, repair kinked lines, or replace the circuit entirely. If not corrected, contact an authorized service provider.

High-Temperature Shutdown

The Travel Air® will shut down if the internal temperature exceeds 132.8°F or 56°C and will not restart until the thermal fuse is replaced. In the event of a device shut down, check that the operating conditions are within the range outlined in this manual. Unplug the unit and replace the thermal fuse (see Maintenance section). If the Travel Air® continues to shut down under normal conditions, contact an authorized service provider.

Chapter 8: Cleaning and Disinfection

 **SAFETY PRECAUTIONS AND IMPORTANT SAFEGUARDS:** To reduce the risk of infection, illness, and injury from contamination or improper use, clean the Phasitron® P5-TA before first use and as soon as possible after each use. Do not submerge the Travel Air®, multimeter, or AC power cord. The Phasitron® P5-TA is for single patient use. Do not share your Phasitron® P5-TA with any other people. Do not mix old and new parts.

 **CAUTION:** Replace the Phasitron® P5-TA breathing circuit according to guidelines provided on the label. Not adhering to time limitations may result in ineffective therapy and Travel Air® malfunction.

 **CAUTION:** No Percussionaire® controller device, single patient Phasitron® P5-TA, or tubing harness should ever be processed in a steam autoclave.

Homecare Cleaning Schedule

Controller	Wipe down weekly or sooner if visibly dirty.
Phasitron® P5-TA	Complete basic cleaning after each treatment. Complete disinfection weekly (after regular basic cleaning of the unit).
Tubing Harness	Complete disinfection and cleaning of the tubing harness weekly.

Travel Air® Controller

Clean and disinfect the surface of the Travel Air® when visibly dirty, weekly, and when transferring between patients (after basic cleaning of the device).

Clean the Travel Air® with a clean, lint-free cloth or paper towel moistened with the cleaner. Use only approved cleaners.

NOTE: Remove the Phasitron® P5-TA tubing harness from the Travel Air® device before cleaning and disinfection.

 **CAUTION:** Do not spray any cleaning solution directly onto the device.

 **CAUTION:** Do not immerse or allow liquids to access the device.


Phasitron® P5-TA Breathing Circuit Kit

NOTE: Phasitron® P5-TA breathing circuit kits are for *SINGLE PATIENT* use only. Do not share your Phasitron® P5-TA with any other people.


⚠ CAUTION: Tubing lines must not be connected during cleaning process.

Disassembly of the Phasitron® P5-TA


Before cleaning, disassemble the Phasitron® P5-TA breathing circuit. Disconnect tubing from the controller device and the Phasitron® P5-TA.


1. Rotate nebulizer bowl towards rear of Phasitron® P5-TA, moving gently past the stop to release.


2. Gently lift nebulizer cap and bowl to separate from Phasitron® P5-TA.


3. Holding the nebulizer cap, twist the nebulizer bowl to remove bowl from cap.
Discard any unused medication in accordance with hospital/institutional protocol and/or medication guidelines.

4. Twist white cap on the rear of Phasitron® P5-TA to remove.

5. Remove cap.

6. Remove sliding venturi and spring from Phasitron® P5-TA body.

Cleaning Procedure

Rinse the Phasitron® P5-TA after each use with sterile water and allow to air dry. Cleaning should be performed daily. Visually inspect interior and exterior of all parts.

1. Before starting the cleaning process, wash hands thoroughly with soap and water or use an alcohol-based hand sanitizer. Visually inspect the exterior of all parts, including tubing, for corrosion, discoloration, pitting, and cracks, as well as cracked or missing O-rings.

 **CAUTION:** Disconnect tubing from Phasitron® P5-TA and the device before cleaning.

 **CAUTION:** Do not wash or submerge the antibacterial filter.

 **CAUTION:** Use of cleaning methods not outlined in these instructions for use may damage the Phasitron® P5-TA and accessories.


2. Use a clean, damp cloth to wipe the exterior of the tubing harness.

3. Thoroughly rinse each of the disassembled parts (except for filter) under warm running tap water for approximately 10 seconds.

4. Use fragrance-free liquid dish soap added to a bowl or basin filled with warm tap water (example: Dawn liquid dish soap in the USA). Do not use "white" liquid dish soaps or antibacterial liquid dish soaps as these may contain additives harmful to parts of the Phasitron® P5-TA kit.

EN

5. Hand wash all parts of the Phasitron® P5-TA and accessories in the warm soapy water. Use a pipe cleaner to clean the monitoring port with warm soapy water.

NOTE: Before cleaning or disinfection, the Phasitron® P5-TA tubing harness must be removed from the Travel Air® device.

6. Use a soft bristle, nylon brush (not supplied) to remove stubborn soil.
7. Rinse thoroughly under running tap water.
8. Inspect prior to reassembly for any remaining soil or damage.

Disinfection Solution Procedure

Use standard household bleach (sodium hypochlorite 5.25%) for disinfection. Mix solution of one-part household bleach and eight parts water. Example: 8 oz bleach to 64 oz water (250 ml bleach to 2 liters water).

⚠ WARNING: Review and follow all safety precautions printed on bleach container.



1. Wear disposable gloves.
2. Pour bleach solution into airtight container large enough to hold 64 oz of water (2 liters).
3. Completely immerse all the Phasitron® P5-TA parts (including mouthpiece and mask) into bleach solution. Do NOT disinfect hydrophobic filter or tubing harness.
4. Let soak for 30 minutes minimum. Do not exceed one hour.
5. Rinse off cold-disinfection solution using sterile or filtered water (less than or equal to 0.2-micron filter). Do not use tap water.
6. Gently shake all parts to remove as much water as possible. Dry with a clean, lint-free cloth.
7. Place on a clean, lint-free cloth and allow to air dry completely.
8. Insert all parts into sealable plastic bag and store in a dry, clean location.

Dry Inside of Tubing







Dry the inside of the tubing with air from the controller device. Connect each tube, one at a time, to the yellow connector on the controller device. Turn on device and run, drying each tube for two minutes, to remove moisture. Hang tubing to air dry.


Phasitron® P5-TA Lubrication After Cleaning and Disinfection

Lubrication is needed only after cleaning and/or disinfection of your Phasitron® P5-TA breathing circuit.

 <p>quick-connect fitting</p>	<p>Wash hands thoroughly with soap and water or use an alcohol-based hand sanitizer before reassembly and lubrication. Lightly coat each quick-connect fitting O-ring with the Percussionaire® lubricant Lubetube, supplied with the Impulsator® device.</p> <p>⚠ CAUTION: Use only Percussionaire®-approved lubricant.</p>
	

Reassemble the Phasitron® P5-TA

	<p>⚠ CAUTION: Do NOT lubricate diaphragm.</p>
	<p>1. Visually inspect O-ring on venturi for cracks or damage.</p> <p>⚠ CAUTION: Do NOT lubricate venturi O-ring.</p>
	<p>2. Visually inspect O-ring on nebulizer cap for cracks or damage.</p> <p>⚠ CAUTION: Do NOT lubricate nebulizer O-ring.</p>
	<p>3. Insert venturi and spring into Phasitron® P5-TA body.</p>
	<p>4. Screw white cap onto rear of the Phasitron® P5-TA body until stop.</p>
	<p>5. Holding the ring of the nebulizer cap, insert and twist nebulizer bowl until it locks.</p>

6. Align the notch, insert the nebulizer housing into the Phasitron® P5-TA body, then rotate to lock.

7. Add mouthpiece, mask, and corrugated tubing to the Phasitron® P5-TA, as needed.

Chapter 9: Maintenance


Home Maintenance Schedule

Service	Time Interval
Pre-use check (see Travel Air® Pre-Use Check section)	When setting up; weekly
Lubrication of connector O-rings and Phasitron® P5-TA	Weekly
Drain water trap	As needed
Replace	Time Interval
Filters	Every 6 months If operating in dirty or dusty environments, replace as needed.
Multimeter batteries	Every 6 months
Phasitron® P5-TA breathing circuit	Every 6 months
Fuses	As needed

Lubrication

After the Phasitron® P5-TA breathing circuit has been cleaned and sanitized, be sure to lubricate the connector O-rings.


Wash hands thoroughly with soap and water or use an alcohol-based hand sanitizer.

1. Coat each connector O-ring fitting with the provided Percussionaire® Lubricant Lubetube.

2. The Phasitron® P5-TA breathing circuit is ready to be reassembled for use. See Cleaning and Disinfection section for reassembly instructions.

Drain Water Trap


Drain water trap by depressing the valve tip at the bottom of the bowl when device is in operation.

Valve tip

Filter Replacement

⚠ CAUTION: The filter elements cannot be washed. Replace filters every six months or sooner if noticeably clogged. Replace filters more frequently when operating in environments with heavy dust or debris.

Water Trap Filter

To replace the Water Trap filter, unscrew the plastic bowl (A) with gasket in a counterclockwise direction by hand and remove the manifold. Unscrew the filter (B) by holding the black portion of the cylinder and turning in a counterclockwise direction to remove from manifold (C). The filter (B) may be cleaned under running water and allowed to air dry.

Water Trap Reassembly

To reassemble water trap, remove handle and place Travel Air® upside down. Replace the filter (B) by gently hand tightening in a clockwise direction into the manifold (C). Hand tighten plastic bowl (A) in a clockwise direction back into manifold.

⚠ CAUTION: Do NOT overtighten; this will damage the gasket. Correct gasket seating will ensure the correct seal is maintained.




NOTE: The Travel Air® water trap filter color may vary. This will not impact the performance of the filter in any way.

Water Trap Gasket

 <p>Gasket</p>	 <p>Gasket placement in water trap Filter color may vary</p>
 <p>Assembled water trap</p>	<p>Always inspect the water-trap gasket when disassembled. If damaged, replace (PRT-B13525-60SL-BL).</p> <p>(See preceding page for water trap reassembly).</p>



Patient Air Filter Replacement

The patient air filter is located behind the water trap. Remove the water trap before replacing the patient air filter. Refer to water trap removal instructions in Chapter 9.


 	<p>1. After removing the water trap, pull out the patient air filter holder.</p>
  	<p>2. Separate the two halves and remove the circular filter.</p>
 	<p>3. Replace the filter (part number PRT-B13608) and insert it into the Travel Air® device. Reinstall the water trap.</p> <p>NOTE: It may be necessary to compress the filter for a proper fit inside the filter holder.</p>

⚠ CAUTION: Do not wash the patient air filter elements. Filters are to be replaced after six months of operation or if noticeably clogged. Heavy dust and debris may cause filters to plug more frequently; maintenance should be provided at increased intervals when operated in such conditions.

Cooling Air Filter

If necessary, trim the replacement blue filter corners for a proper fit (Percussionaire® part number PRT-B13091).


To replace the cooling air filter behind the grill, open the grill by turning the latch counterclockwise using a coin or screwdriver.


EN

Fuses

 **WARNING:** The power cord must be disconnected when changing the fuses.


Switch Fuse

The switch fuses can be checked and/or changed by pulling the holder straight out.

Thermal Fuse


The thermal fuse can be checked/changed by twisting and pulling the holder straight out.

 **CAUTION:** If the fuses blow immediately upon restarting, contact an authorized service provider.

PDM (Multimeter) Batteries

Side view of PDM

Install batteries into PDM.

• USB serial port only used for firmware upgrades

 A Low Battery indicator is displayed when battery capacity is nearing depletion.

Changing PDM Batteries

1. Press on the PDM's bezel and twist counterclockwise approximately 20 degrees.
2. Gently pull on the multimeter to remove it from the housing.
3. Remove the two old batteries.
4. Install two new batteries. Note that the positive terminals face the same direction. Wait 30 seconds until the screen turns off.
5. Install the PDM back into the housing and twist clockwise until the stop is felt.
6. Refer to System Description section: POST mode to verify display operation.

NOTE: Do not install PDM until the POST check is complete, and the screen is blank, indicating Sleep mode.

EN

NOTE: The PDM has a USB serial port that is used for manufacturing, calibration, and firmware upload. It is not enabled during normal operation. It is not intended to be used by the patient or in normal operation of the PDM.

Chapter 10: Accessories

Use only accessories recommended by Percussionaire® for maximum performance and safety.

Description	Order Number
Phasitron® P5-TA kit	P5-TA-5
Mouthpiece (10 pack)	PRT-B20018
Battery CR123A (2 pack)	PRT-B13350
Lubrication Lubetube	A50060
A55122 Travel Air® Accessory Kit *Contents not sold individually	Kit Includes:
Cooling Air Filter	B13091
Patient Air Filter	B13608
Fuse, Compressor	B13326-3
Fuse, Power Switch	B13327
Gasket, Water Trap Bowl	B13525-60SL-BL
Plug, Cap	B13536

Chapter 11: Troubleshooting

Problem	Possible Cause	Remedial Action
The Travel Air® will not start.	The device is not connected to a power source. A fuse is blown.	Plug the device into power source. Replace the fuse. If the fuse blows immediately upon starting the device, contact an authorized service provider.
The Travel Air® fails to maintain pressure.	There is a leak at the patient mask or airway. There is an internal leak.	Adjust the patient mask or airway to stop the leak. Contact an authorized service provider.
The Phasitron® P5-TA will not percuss.	The Phasitron® P5-TA was not assembled correctly.	Refer to the assembly section of the Travel Air® User Manual and follow Phasitron® P5-TA assembly instructions.
The Phasitron® P5-TA nebulizer will not aerosolize.	The nebulizer baffle is loose (small brown tube inside bowl). There is no gas flow from the yellow line.	Inspect and reseal nebulizer baffle on stem. Contact an authorized service provider.
The pressure and/or frequency of percussions is not correct.	The Travel Air® is not calibrated correctly.	Contact an authorized service provider.
The multimeter is inoperative.	The red line is loose or leaking. The red line filter is blocked.	Make sure tubing connections to Phasitron® P5-TA are correctly attached. Clear or replace the red line filter.
The Travel Air® shuts down.	There is an obstruction in the tubing or within the Travel Air® device.	Check the tubing for obstructions. Contact an authorized service provider.

NOTE: For issues that are not covered above, please contact Percussionaire® or authorized service center.

Chapter 12: Service and Repair

Percussionaire® recommends an annual Preventive Maintenance (PM) for each device. An annual PM consists of a thorough cleaning, functional evaluation, and, if necessary, recalibration. Intervention by an unauthorized individual or repair maintenance facility will cause the immediate expiration of the clinical readiness of the device. Return the device to Percussionaire® or an authorized factory service center for repair or annual preventive maintenance.

Disposal of Equipment

At the end of useful life of a unit, disposal should be in accordance with local, state, federal and international laws. The battery is considered electronic waste and must be disposed of according to local regulations. Follow local governing ordinances and recycling plans regarding disposal or recycling of the battery.

Returns

To return the Travel Air® device, contact your distributor. To pack the device, use the black carry bag with the enclosed foam. This will protect the water trap during shipping.

Chapter 13: Technical Specifications

Travel Air® Controller

Settings	Specifications
Aerosol flow	Up to 25 LPM
Percussion frequency	60-500 pulses per minute, 1-8.3 Hz
Pulse/interval ratio	Automatic
Operating Conditions	Specifications
Operating range	Temp., 5°C to 40°C (41°F to 104°F) Humidity 5%-95%
Storage and transport range	Temp., -20°C to 60°C (-4°F to 140°F) Humidity < 93% non-condensing
Power input	100-240 VAC, 50/60 Hz, 3 amp
Run Time	20 minutes
Display/Output	Specifications
Airway pressure	Digital display, 0-99 cmH ₂ O, accurate to +/-1 cmH ₂ O/hPa
Pulse Frequency	Digital display, 50-999 pulses per minute
Pulse Amplitude	Digital display, 1 to 99 cmH ₂ O/hPa, accurate to +/-1 cmH ₂ O/hPa
Treatment time duration	Digital display
Session usage time	Digital display, maximum 59:59 (minutes:seconds)
Alarm Condition	Indication
Airway pressure over 90-100 cmH ₂ O	Audible indicator
Compressor overload	Audible indicator/Shut down
Compressor over temperature	Shut down
Compressor	Specifications
Pressure	40 PSI
Flow	30 LPM
Filters	Description
Patient air filter	Compressor air intake filter
Cooling air filter	80 mm intake filter
Dimensions/Weight	Specifications
Dimensions (W x H x D)	26.67 cm x 25.4 cm x 25.4 cm (10.5" x 10" x 10")
Weight	8.98 kg (19.8 lbs)
Required Maintenance	Filter changes every 6 months, or sooner if required
Service Life	6 years from date of manufacture

Phasitron® P5-TA Specifications

Size	13.5 mm x 17 mm (5 ¼" x 6 ¾")
Weight	123 g (0.27 lb)
Operating Range	Temp., 0°C to 49°C (32°F to 120°F) Relative humidity range, 15% to < 90% non-condensing
Storage and Transport	Temp., -40°C to 60°C (-40°F to 140°F)
Rate Range	0-999 pulses per minute
Pressure Range	0-150 cmH ₂ O/hPa
Aerosol Flow	.75 cc per minute
Safety Valve Release	30-50 cmH ₂ O/hPa
Red Line Filter	1-3 micron hydrophobic
Disposal	Recycle according to local laws
Service Life	6 months or 540 uses, whichever is less
Shelf Life	2 years from date of manufacture

Digital Multimeter (PDM) Specifications

Size	73 mm diameter (2.87-inch diameter)
Mass	165 g (0.36 lb)
Storage and Transport	Temp., -20°C to 60°C (-4°F to 140°F) Humidity < 93% non-condensing
Display	128 x 64-pixel FSTN chip-on-glass LCD with reflector
Fault Detection	Independent hardware and software watchdogs
Serial Port	USB (Manufacturing, firmware upgrade use only)
Thermal	-40°C to +60°C (-40°F to 140°F)
Rate Range	50-999 pulses per minute
Pressure Range	1-150 cmH ₂ O/hPa
Pressure Resolution	1 cmH ₂ O/hPa
Pressure Accuracy	Greater of ± 0.5% of reading or 1 cmH ₂ O/hPa
Battery Type	CR123A 3.0V (2)
Battery Duration	3,250 operational hours at 35°C (95°F)
Shelf Life	3.5 years at 35°C (95°F)

EN

EMC Warnings

Medical Electrical Equipment needs special precautions regarding electromagnetic compatibility (EMC), and needs to be put into service according to the EMC information provided in this manual. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

The Travel Air® should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, run the Travel Air® before treatment begins to verify the device operates normally in the desired configuration.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The Travel Air® is intended for use in the electromagnetic environment specified below. The user of the Travel Air® should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purpose.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

The Travel Air® is intended for use in the electromagnetic environment specified below. The use of the Travel Air® should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±8 kV air	±8 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/burst IEC 61000-4-4	±2 kV for supply mains	±2 kV for supply mains	Mains power quality should be that of a typical home or hospital environment.
Electrical surge immunity IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	± kV line to line ± kV line to earth	Mains power quality should be that of a typical home or hospital environment.
Voltage dips, short interruptions and voltage on power supply input lines IEC 61000-4-11	Mains Voltage 30% for 25 cycles 60% for 60 cycles 95% for 0.5 cycles 95% for 5 sec	Mains Voltage 30% for 25 cycles 60% for 60 cycles 95% for 0.5 cycles 95% for 5 sec	Mains power quality should be that of a typical home or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device to be powered from an uninterrupted power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical home or hospital environment.

NOTE: *UT is the s.c. mains voltage prior to application of the test level.*

The Travel Air® is intended for use in the electromagnetic environment specified below. The user of the Travel Air® should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 610004-6	3V RMS 80 MHz to 2500 MHz	3V RMS 80 MHz to 2500 MHz	Mains power quality should be that of a typical home or hospital environment. Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.7 \cdot \sqrt{P}$ $d = 1.7 \cdot \sqrt{P}$ $d = 1.2 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \cdot \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2500 MHz	3V/m 80 MHz to 2500 MHz	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended Separation Distances

It is important to follow the guidelines for recommended separation distance between portable and mobile RF communications equipment and Travel Air®.

The Travel Air® is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Travel Air® can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Travel Air® as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Power Output of Transmitter (Watts)	Separation Distance According to Frequency of Transmitter (Meters)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	0.37	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

Note 3: An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range of 80 MHz and 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

Note 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


EN

Chapter 14: Limited Warranty

Percussionaire® warrants that the Travel Air® shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of one year from the date of first use (proof of delivery will be required). If the product fails to perform in accordance with the product specifications, Percussionaire® will repair or replace – at its option – the defective material or part. Percussionaire® will pay customary freight charges to and from Percussionaire or an authorized Percussionaire® service center. This warranty does not cover damage caused by non-approved cleaning or sterilization, accident, misuse, abuse, alteration, and other defects not related to material or workmanship. Percussionaire® disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product.

NOTE: In the event that the operation of a Percussionaire® product is in any way adversely affected by using components other than those designed, manufactured, or approved by Percussionaire® Corporation, Percussionaire® shall not be liable under this warranty with respect to such product.

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