ThevoautoActiv

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Dear Customer,

Thank you for choosing our product and for the trust you have placed in us. In order for you and the user to enjoy your new product safely, conveniently, and comfortably, we kindly ask that you read the Instructions for Use. Doing this will also help avoid damage, injury, and incorrect operation and handling of the product. Please ensure that a copy of the Instructions for Use remains with the product so that every user has access to it.

If you have any questions, please contact your authorised distributor or us directly.

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We reserve the right to make technical changes. This can lead to differences between your product and how it is presented in brochures and other media.

Table of Contents

Figures	3
Preface / Manufacturer	18
Labels / Model Identification	19
General Remarks / Safety Instructions	20
Intended Purpose	20
– Intended Use / Application Areas	20
– Indication / Patient Target Group / Intended Users	20
– Risks Associated with Use and Contraindications	21
Product Description	21
– Scope of Delivery / Standard Equipment	21
– Options	21
– Putting into Service	22
– Teaching the User	22
Adjustment Options / Features	22
– Accessories	23
– Technical Data	24
– Material Specifications	25
Cleaning, Maintenance, and Repair	25
– Cleaning, Care Instructions, and Disinfection on a Daily Basis	25
– Maintenance and inspection by the distributor / Maintenance Schedule	26
– Spare Parts	26
– Storage / Disposal	26
EU Declaration of Conformity	26
Warranty	27
Period of Use / Service Life	27
Reuse	27

Labels / Model Identification

Label	Meaning																								
<p> Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG Walkmühlenstraße 1 - 27432 Bremervörde - Germany ThevoautoActiv 135 90x200 Inko Thevo 64851/0100 9800648510100 MD CE 2021-06-03 UDI-DI (01)04034089082813 UDI-PI (21)1901234 SN 1901234 </p>	<p>Barcode Label</p> <p>Position on the Product (Fig. 1 + 2)</p> <ul style="list-style-type: none"> - Manufacturer - Product Name - Product Group - Item Number short / long - MD = Medical Device - Manufacture Date - Serial Number 																								
<div style="display: flex; justify-content: space-around;"> <div style="width: 45%;"> <p>Thomashilfen ThevoautoActiv Kopfende hier / head end here</p> <table border="1"> <thead> <tr> <th>Service</th> <th></th> </tr> </thead> <tbody> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </tbody> </table> <p>www.thomashilfen.de</p> </div> <div style="width: 45%;"> <p>Thomashilfen ThevoautoActiv < 135 kg</p> <table border="1"> <thead> <tr> <th>Service</th> <th></th> </tr> </thead> <tbody> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </tbody> </table> <p>www.thomashilfen.de</p> </div> </div> <p> ThevoautoActiv 135 90x200 Inko Thevo 64851/0100 9800648510100 MD CE 2021-06-03 UDI-DI (01)04034089082813 UDI-PI (21)1901234 SN 1901234 </p>	Service												Service												<p>Maintenance Label + Barcode Label</p> <p>Position on the Product (Fig. 1 + 2)</p> <p>The maintenance label documents who did the regularly scheduled maintenance work and when.</p> <ul style="list-style-type: none"> - Mattress (Fig. 1) - Suspension framework (Fig. 2)
Service																									
Service																									

Remarks

General Remarks

- Read the Instructions for Use carefully before first use. If you have difficulty reading, please have them read to you. If the copy of the Instructions for Use is missing, you can download one from www.thomashilfen.de. The contents of the digital Instructions for Use can be enlarged.
- Check the packaging for any transport damage and to be sure the product is in flawless condition.
- Changes to the product are only permitted within the scope of the intended adjustment options.
- Cleaning and daily care can be provided by relatives and caregivers. Maintenance and inspection are done by the distributor.
- All labels have to remain on the product and should not be tampered with.
Exception: The maintenance label has to be updated and completed after each maintenance.

Safety Instructions

- Please do not overload the product and note the maximum user weight and the maximum load capacity (see technical data).
- The product should only be used by one user at a time. No other individuals should be placed in or on the product.
- Heavy perspiration, high humidity, and extreme changes in temperature can have an impact on the foam in the mattress. Pay attention to changes such as sagging in the mattress.
- If there is a serious and unexpected accident with the product, please inform our company immediately and the responsible authorities.

Intended Purpose

Intended Use / Application Areas

The ThevoautoActiv micro-stimulation system is a class 1 medical device.

The system is suitable for use in home care and in professional care settings. It consists of a mattress and a system of suspensions with the MiS micro-stimulation effect that promotes natural movement and the user's body awareness. The system allows users to lie a specific way; it can be complemented as required by other positioning aids. Adapters can target and relieve pressure from vulnerable areas of the body. Special air channels in the foam core allow for very good air circulation so that moisture can be removed. The micro-stimulation system is intended to treat and prevent pressure sores and to reduce pain caused by lying down.

The system replaces a standard mattress and goes into existing bed frames with the reclining surface and is also used in beds in care settings. The adjustment options of the existing bed frame are not limited by the micro-stimulation system. The system does not replace manually changing positions, but the system can lengthen the lying intervals. It is suitable for use in combination with digital assistance systems in care settings.

Any use that goes above and beyond this does not comply with the intended use.

Indication / Patient Target Group / Intended Users

The micro-stimulation system is suitable for adolescents and adults with limited mobility and / or impaired body awareness, circumstances that adversely impact changes or shifts to body position. Furthermore, the product is used if the user has or is at risk of having pressure sores on the skin and subcutaneous tissue, such as on the back, buttocks, or heels.

This can definitely be used if there are illness- or age-related difficulties falling asleep, sleep disturbances, and impaired body awareness during sleep.

Risks Associated with Use and Contraindications

If used properly, risks can be eliminated during use. There are no known contraindications.

Product Description

Scope of Delivery

Your product will be delivered in several boxes and must be assembled when commissioned. A copy of the Instructions for Use is enclosed with the product.

Standard Equipment

Mattress core; suspension framework; cover (jersey or incontinence)

Components

Active side element (Fig. 3)

The active side element does not only connect the single flexible fibreglass lath, but also offers the possibility to integrate ThevoautoActiv into common nursing and hospital beds. You are able to carry on using the head and foot adjustment. The system is simply laid on the existing slatted frame or inlay frame. Due to the flexible construction you can still use the head and foot adjustment as well as the height adjustment of the bed.

Adaptor (Fig. 4)

The adaptors are situated in the active side element (ASE). Through these adaptors the suspension travel of the flexible fibreglass lath can be restricted as required. Thus, through various adjustments body parts at risk can be alleviated additionally. By taking the adaptors completely off, the highest possible alleviation is reached.

Flexible fibreglass lath (Fig. 5)

These laths are made from a durable material by innovative glass fibre technology (pultrusion). Beneath a long durability they ensure an optimum lying comfort, which has a favourable effect on the convalescence of the patient.

Wing suspension (Fig. 6)

3 or 4 wing suspensions situated on one flexible fibreglass lath optimally adapt to the body contours of the patient. An anatomically correct body alignment is ensured and sleep is influenced positively. Each wing suspension can be moved or removed in order to support the individual preventive and therapeutic treatment in an optimum way.

Options

Upon request the product can be made in the following sizes:

- Width options from 70 cm to 140 cm (27.6" to 55.1")
- Length options from 180 cm to 220 cm (70.9" to 86.6")

The mattress is available for the following weight categories:

- **ThevoautoActiv 50** – from 30 kg to max. 50 kg patient weight (66.1 to 110.2 lbs)
- **ThevoautoActiv 100** – from 50 kg to max. 100 kg patient weight (110.2 to 220.5 lbs)
- **ThevoautoActiv 135** – from 100 kg to max. 135 kg patient weight (220.5 to 297.6 lbs)

If the patient is in the upper end of a weight category (+/- 2 kg / 4.4 lbs), the next higher weight category should perhaps be chosen.

If the patient's body has an atypical distribution of weight, the next weight category should perhaps be chosen.

Putting into Service

- The product can be used in any bed. Put the mattress with the red side up onto the slatted frame or your bed's supporting structure.
- The mattress can also be used without any slatted frame. In this case, there has to be a support surface of at least 13 cm / 5.1" in width on the long sides.
- Please note: For beds with side rails, side rail risers might be need.

Assembly

Please note, that the ThevoautoActiv mattress will be delivered in an extra carton. The system will be delivered in two parts (Fig. 7).

1. Please put the two parts close together, so that the side elements fit into each other. You should always pay attention on the direction marking (Fig. 8).
2. The connection thorns have to be placed tight and flush into the side elements at the connection points 1 to 6 (Fig. 9 + 10).
3. Lay the under frame on the existing slatted frame or inlay frame of your bed. Please pay attention to the label for the head and foot section.
4. Now lay the mattress on the under frame.
5. The black fixation flaps of the mattress cover in the head and foot section must be flapped around the under frame and closed with the zipper (Fig. 11).
6. The under frame of ThevoautoActiv is safely fixed to the bed frame with the 6 belts. Thread the belts through the bottom holes of the active side element and then around the tubular frame of the bed. Thread the belt through the buckle and pull it tight. To release the belt, open the buckle and pull out the belt (Fig. 12).
Attention: If you have a height adjustable bed, please fix the belts to the lifting system.
7. The system can now be individually adjusted to meet the needs of the respective patient.
8. Check the settings after 24 hours at the latest and carry out eventually necessary changes.

Teaching the User

The user is taught how to use this device by trained staff who will reference these Instructions for Use.

Adjustment Options / Features

Standard setting at delivery

At the delivery of ThevoAutoactiv the adaptors and flexible fibreglass laths are set as shown on Fig. 13.

Adjusting the adaptors

In order to adjust the adaptors you have to grab under the mattress and pull the flexible fibreglass-track slightly upwards. Now you can turn the adaptor to the required direction or rather remove it completely. Please ensure that the adaptor grips tight into the stabiliser after every adjustment. Only in this case an unintended turning can be avoided (Fig. 14).

Adjusting the flexible fibreglass lath

First you will have to take the adaptor out off the side element. Push the flexible fibreglass lath out off the stabiliser (the stabiliser is situated under the flexible fibreglass lath in the side element) (Fig. 15). Now you can take the FGF- track out off the side element. Put the stabiliser in the desired position of the side element and slide in the flexible fibreglass lath (Fig. 16). Put in the adaptor again and turn it into the required position.

Advanced improvements

We will show you nine setting options of which you shall choose the right one for your patient in order to accordingly set ThevoautoActiv. These recommendations are based on the results of our carefully checked experiences and for that reason they shall be viewed as average values. Please keep in mind that we can only give recommendations for the individual. The individual adjustment should be made because of the clinical picture. We recommend a regular check of the settings to ensure also the correct treatment in the progressing therapy. When adjusting the adaptors please act step-by-step from the head to the foot and compare your personal setting with the corresponding picture in this Instructions for Use. The adjustment of the adaptors should be carried out as described on.

Important hint: Only remove the adaptors and not the flexible fibreglass laths.

Setting recommendation for a body height of 150 - 169 cm (59.1 - 66.5")

- Light weight, light physique (Fig. 17)
- Medium weight, normal physique (Fig. 18)
- Heavy weight, strong physique (Fig. 19)

Setting recommendation for a body height of 170 - 184 cm (66.9 - 72.4")

- Light weight, light physique (Fig. 20)
- Medium weight, normal physique (Fig. 21)
- Heavy weight, strong physique (Fig. 22)

Setting recommendation for a body height of 185 - 200 cm (72.8 - 78.8")

- Light weight, light physique (Fig. 23)
- Medium weight, normal physique (Fig. 24)
- Heavy weight, strong physique (Fig. 25)

Important hint: If the person is lighter than 30 kg (68 lbs) please get in contact with us.

Accessories

Extra mattress cover (Fig. 26)

An extra mattress cover made of jersey or incontinence material is available if the original mattress cover has to be cleaned or discarded.

ThevoFloor – Floor bedding system (Fig. 27)

ThevoFloor is suitable for use at home and in professional care settings and allows patients to lie low to the floor. This reduces the need for measures which infringe on the freedom of the user and prevents falls out of bed. The system is comprised of:

- two wall parts (foam with incontinence cover)
- connecting belt
- cover strip for the hook and loop fastener
- A mattress and a system of suspensions with the MiS micro-stimulation effect are used that promote natural movement and the user's body awareness.

The system can be designed to suit the user. A standard Thevo mattress measuring 140 x 200 cm / 55.1 x 78.7" fits. Other sizes are available on request.

ThevoFloor is used in the corners of rooms. The Thevo mattress selected has to be appropriate for the indication and weight of the patient in question.

Put both wall parts on the walls in the corner of the room. Make sure that the hook and loop strip on the front of the short, top-wall part faces the hook and loop strips on the front of the longer side wall part (Fig. 28). Connect the wall parts firmly with the hook and loop strip. The cover strip can be put on the hook and loop fastener.

Pull the supplied connecting belt through the loops in the cover of the wall parts (Fig. 29).

Put the mattress red side up in the corner of the room. The top of the mattress should face the short part of the wall. Bring the connecting belt around the mattress. Adjust the length of the belt so that the wall parts are securely connected with the mattress (Fig. 30).

ThevoFloor wall parts	Top part	Side part
Length	140 cm / 55.1"	208 cm / 81.9"
Height	80 cm / 31.5"	80 cm / 31.5"
Width	8 cm / 3.1"	8 cm / 3.1"

Rescue belt system (Fig. 31)

With the rescue belt system bedridden patients can be removed quickly from danger. The rescue belt system is integrated into the cover of the mattress. The patient is secured by two hook and loop safety belts (at the chest and thigh). In an emergency, the hook and loop belts are pulled over the mattress, bedding, and patient and are connected. This way the patient and the bedding are secured quickly and easily to the mattress. The patient's arms are under the covers so that they are not injured during transport.

Using the pull straps at the top and foot of the bed the mattress, patient, and bedding are pulled out of the bed at a 90° angle. This way the patient can be evacuated horizontally and vertically from danger. The mattress slides very easily on all floor types. There are two lifting straps on each long side which allow emergency personnel to take the patient without transferring them from the bed.

The fabric of the belts is 100% polyester and is free of chemicals. This makes them non-irritant and also suitable for patients with allergies.

Maximum load capacity: up to 150 kg / 330.7 lbs

All users have to be familiar with how the product works. This way in an emergency the user can be removed quickly from the product.

Technical Data

ThevoautoActiv 50

Overall dimensions (w x l)	90 x 200 cm / 35.4 x 78.7"	100 x 200 cm / 39.4 x 78.7"
Height	19 cm / 7.5"	19 cm / 7.5"
Weight	21,4 kg / 47.2 lbs	23,2 kg / 51.1 lbs
max. user weight	50 kg / 110.2 lbs	50 kg / 110.2 lbs

ThevoautoActiv 100

Overall dimensions (w x l)	90 x 200 cm / 35.4 x 78.7"	100 x 200 cm / 39.4 x 78.7"
Height	19 cm / 7.5"	19 cm / 7.5"
Weight	22,4 kg / 49.4 lbs	24,3 kg / 53.6 lbs
max. user weight	100 kg / 220.5 lbs	100 kg / 220.5 lbs

ThevoautoActiv 135

Overall dimensions (w x l)	90 x 200 cm / 35.4 x 78.7"	100 x 200 cm / 39.4 x 78.7"
Height	19 cm / 7.5"	19 cm / 7.5"
Weight	24,2 kg / 53.4 lbs	26,3 kg / 58 lbs
max. user weight	135 kg / 297.6 lbs	135 kg / 297.6 lbs

Material Specifications


The open-pored and flame retardant cold foam provides moisture absorption and release and thermal regulation making the foam highly durable. The covers are flame retardant.

Cleaning, Maintenance, and Repair**Cleaning, Care Instructions, and Disinfection on a Daily Basis****Foam Mattress**

The mattress can be autoclaved making it suitable for reuse.

Jersey Cover

The jersey cover is breathable and can be easily unzipped and removed for washing. Please also note the following care instructions.

**Incontinence Cover**

The incontinence cover can be cleaned by wipe-down disinfection. The incontinence cover can be easily unzipped and removed for washing. Please also note the following care instructions.

**Protective Cover for the Suspension framework**

The black cover for the suspension framework can be easily unzipped and removed for washing. Please also note the following care instructions.

**Suspension framework**

Simple Cleaning: The bed base can be damp wiped with non-corrosive household cleaning agents or washed with a wet cloth. Make sure the bed base can dry well after. Cleaning with an alcohol-based disinfectant is also possible.

Disinfection

The mattress and suspension framework can be cleaned by steam vacuuming. This has to be done to disinfect and before each reuse.

We recommend two disinfecting methods recognised by the Robert Koch Institute:

- the fractionated vacuum method (vacuum-steam vacuum method) up to a maximum temperature of 110° C.
- the steam flow method up to a temperature of 110° C.

Important note: Thermal treatment above 110° C leads to premature ageing and deforming of the plastic parts which is why it shouldn't be done. When using chemical disinfectants always follow the information from the manufacturer to ensure they are suitable for plastics.

Maintenance and inspection by the distributor

Preparation

A copy of the Instructions for Use has to be enclosed with the product. If one is not available, request a copy from the distributor. Familiarise yourself with how the product works. If you are not familiar with the product, read the Instructions for Use before testing it. Clean the product before starting to test it. Please note the care instructions in the Instructions for Use and product-specific test instructions. If you have further questions, please contact the staff member you dealt with or Thomashilfen directly.

Service and repairs to the product should only be done by qualified staff. It is the user’s responsibility to ensure that the prescribed inspections / maintenance are done regularly and in a timely fashion. If damage is apparent, the user has to actively inform and have the damage fixed immediately by qualified staff.

Maintenance Schedule – Annual Inspection

	Area	Function	Damage	Replace
1	Cover	👁 / 🔄	👁 / 🔄	
2	Foam core	👁 / 🔄	👁 / 🔄	
3	Suspension framework	👁 / 🔄	👁 / 🔄	
4	Wing suspensions	👁 / 🔄	👁 / 🔄	
5	Suspension blocks	👁 / 🔄	👁 / 🔄	
6	Labels		👁 / 🔄	
done by:			on:	
👁 = Inspection 🔄 = Reuse ✓ = Test completed				

Spare Parts

Only original spare parts should be used for repairs. A list of the spare parts is available at www.thomashilfen.de. To receive the right spare parts or for possible technical recalls, we normally require the serial number / UDI-PI and the exact UDI-DI.

Removal / installation of spare parts should only be done by qualified staff! In the event of returns, please make sure the product is securely packaged for transportation.

Storage

Store the product in a safe place when it is not being used. Do not lay or put any heavy objects on the product while it is being stored. Protect the product from weather like rain, intense sunlight, and frost. If the product is wet, wipe the metal and plastic parts dry and let the padding dry before storing the product.

Disposal

If the product is no longer usable due to its condition or it has reached the end of its service life, it can be recycled at the local waste disposal company. Please refer to the local guidelines for recycling and disposal.

EU Declaration of Conformity

We hereby declare as the party solely responsible that this product and any versions and accessories comply with the basic regulations. The current EU Declaration of Conformity can be downloaded at www.thomashilfen.de or requested from us directly.

Warranty

Thomashilfen offers you the two-year warranty as required by law for this product. The warranty period starts when the item is handed over. The warranty includes all claims which impact function. Not included is damage which is the result of improper use (e.g. overloading) and natural wear and tear.

Period of Use

The product has a period of use of up to 10 years if used correctly and serviced regularly and professionally as per the information from the manufacturer. The product can be used beyond this point if it is in safe and technically sound condition. Time spent in storage at the distributor does not count towards the period of use.

The theoretical period of use does not equate to a guaranteed service life. It depends on the intensity of use and always undergoes an individual assessment by the distributor.

The period of use does not apply to wear parts like covers, the mattress core and so on.

Service Life

The product has a service life of 10 years or 2 reuses (3 users) if the requirements stipulated for the period of use are met. The service life consists of the period of use and the storage times between reuses.

Reuse

The product is suitable to be reused. Before being passed on to the next user the product always has to be carefully cleaned and disinfected. The condition has to be checked by trained staff for any wear and damage. Any issues have to be remedied to ensure the continued safe use of the product. A copy of the Instructions for Use is a component of the product and has to be passed on to the new user.

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