INSTRUCTIONS FOR USE

TheraKair Visio

Mattress Replacment System

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Benutzerhandbuch · Handleiding · Manuel d'utilisation · Manuale d'uso · Manual del Usuario ·
Brugervejledning · Bruksanvisning

Matratzensersatzsystem · Matrasvervangend ligsysteem · Système de remplacement de matelas
Sistema sostitutivo del materasso · Sustitutivo de Descanso (SSD) · Madrasudskiftningsystem
Madrassersättningsystem

User Manual 3

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Handleiding..... 83

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Manuale d'uso 163

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Brugervejledning..... 243

Bruksanvisning 283

WARNING

To avoid injury, always read this Instructions for Use and accompanied documents before using the product.

Mandatory to read the Instructions for Use.

Design Policy and Copyright

® and ™ are trademarks belonging to the Arjo group of companies.

© Arjo 2019.

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. The content of this publication may not be copied either whole or in part without the consent of Arjo.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

Arjo hereby disclaims all express or implied warranties, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Arjo product(s) described in this publication. Under no circumstances shall Arjo be liable for any indirect, incidental, or consequential damages and expenses, including damages or injury to person or property, due in whole or in part to the use of the product other than those for which disclaimer of warranty or limitation of liability is expressly prohibited by specific, applicable law. No person has the authority to bind Arjo to any representation or warranty except as specifically set forth in this paragraph. Descriptions or specifications in Arjo printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Information in this publication may be subject to change at any time. Consult the manufacturer for updates.

Important Information for Users

In order for Arjo products to perform properly, Arjo recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by Arjo. Contact Arjo for information regarding maintenance and repair.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards. To avoid the risk of electric shock, this product must be connected to a grounded power receptacle.

Warnings

- **WARNING:** Use or non-use of restraints, including side rails, can be critical to patient safety. Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. See related Safety Information.
- Fluids remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and staff. Ensure that there is no moisture in or near the power inlet, power switch and power plug before reconnecting the power supply.
- Do not use sharp instruments to open packaging. Damage could result to mattress, cover sheet, therapy unit or mattress hose.
- The control screen should only be operated by finger; using pens or pointing devices will damage the screen and affect future use.

Specific indications, contraindications, warnings, precautions and safety information exist for Arjo's therapeutic support surface products. It is important for users to read and familiarize themselves with these instructions and to consult the treating physician prior to patient placement and product use. Individual patient conditions may vary.

Notice

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltage.

TABLE OF CONTENTS

Disclaimer of Warranty and Limitation of Liability	4
Important Information for Users.....	4
Warnings	4
Notice	4
Introduction.....	8
About the TheraKair Visio Mattress Replacement System	8
System Features	8
Low Air Loss.....	8
Pulsation Therapy	8
Scheduled Pulsation	8
Pressure Relief Therapy	9
Reduced Heel Area Pressure.....	9
InstaFlate™ Function	9
Touch Screen Interface	9
Individualized Patient Support.....	9
Quick Release Air Supply Connector.....	9
Warmer	9
Reusability.....	9
Indications	10
Contraindications.....	10
Intended Care Setting	10
Connecting the System to Other Devices	10
Devices That May Be Attached to This System	10
Devices to Which This System May Be Attached	10
Risks and Precautions.....	10
Mattress Replacement	10
Transfer	10
Side Rails and Restraints.....	10
Patient Migration	11
Use With Other Devices.....	11
Protection Against Hazards.....	11
Fluids.....	11
Power Cable.....	11
Safety Information	11
Patient Size and Weight.....	11
Air Intake	12
Side Rails / Patient Restraints	12
Bed Frame	12
Bed Height	12
Brakes	12
Head of Bed Elevation	13

- Patient Entrance / Exit 13
- Skin Care 13
- No Smoking in Bed 13
- Therapy Unit..... 13
- Interference 14
- General Protocols 14
- End of Life Disposal 14
- Customer Contact Information 14
- Preparation for Use 15
 - Check Out Procedures 15
 - Preparing the System for Use 15
 - Fitting the Mattress to the Bed 15
 - Fitting the Mattress Cover Sheet..... 16
 - Preparing the Therapy Unit 17
- Operation..... 18
 - Power-Up Procedures 18
 - Initial Power-Up..... 18
 - Power-Up After an Improper Shutdown 18
 - Power-Down Procedure 18
 - Performing CPR (Cardiopulmonary Resuscitation)..... 19
 - Control Screen Display..... 19
 - Control Screens..... 20
 - Home Screen 20
 - Utilities Screen 21
 - About Screen 21
 - Settings Screen..... 22
 - Alarms Screen..... 22
 - Language Screen..... 23
 - Access Code Screen..... 23
 - Pulsation Screen 24
 - Warming Screen 24
 - InstaFlate* Screen 25
 - Firmness Screen..... 25
 - Presets Screen..... 26
 - Quick Reference Guide 26
- Alarms and Alerts 27
 - System Integrity Alarm 27
 - Air Pressure Alarm 27
 - Over Temperature Alarm 28
 - Therapy Inactive Alarm 28

- Mattress Not Connected Alert28
- Patient Placement29
 - Preparing for Patient Placement29
 - Patient Transfer to *TheraKair Visio* MRS29
 - Patient Specific Adjustments30
 - Foot End Interface.....30
 - Pulsation Adjustment.....30
 - Warming Adjustment.....31
 - Firmness Adjustment.....31
 - Patient Transfer from *TheraKair Visio* MRS31
 - Patient Transport on the *TheraKair Visio* MRS31
- Nursing Care32
 - Bathing the Patient.....32
 - After Bathing the Patient32
 - Skin Care.....32
 - Incontinence / Drainage32
 - Mattress Firmness Adjustment33
 - Bedpan Replacement.....33
 - Bedpan Removal.....33
- Care and Cleaning.....34
 - Daily Care and Cleaning of *TheraKair Visio* MRS While in Use.....34
 - Mattress34
 - Therapy Unit Screen34
 - Weekly Care and Cleaning of *TheraKair Visio* MRS While in Use.....35
 - TheraKair Visio* MRS Cover Sheet.....35
 - Inflatable Air Cushions35
 - Power Cable.....35
 - Periodic Maintenance of the Therapy Unit36
 - Safe Transport and Storage of the *TheraKair Visio* MRS.....36
 - Transport and Storage of Individual System Parts.....36
 - Transport of System While Fitted to Hospital Bed36
- Specifications37
 - Classification37
 - Mattress Overlay37
 - Therapy Unit.....37
 - Electrical Data37
 - Transport and Storage Conditions38
 - Operating Conditions38
- Replacement Part Numbers38
- Electromagnetic Compatibility (EMC).....39
- Definition of Symbols Used42

INTRODUCTION

About the TheraKair Visio Mattress Replacement System

The TheraKair Visio™ Mattress Replacement System uniquely combines low air loss with a nylon fabric with moisture vapor permeable backing and pulsation therapy to help provide the ideal environment for the skin of the patient.

The *TheraKair Visio* MRS can reduce pressure, shear and friction through low air loss and through the use of nylon fabric with moisture vapor permeable backing, enables skin moisture and temperature to be managed. Pulsation therapy mimics the body's natural movements, encouraging lymph and blood flow for increased oxygenation to the skin which may help reduce the risk of skin damage.

The system features an intuitive touch screen for easy operation and has been designed to be silent for patient comfort.

System Features

The *TheraKair Visio* MRS offers the following features and advantages over conventional air filled overlays or hospital bed mattresses:

Low Air Loss

- May help reduce pressure in comparison to leading alternating pressure mattresses
- Conforms to bony prominences, minimizing deformation of body tissue
- Allows the patient to sink into the mattress, evenly distributing pressure
- Helps reduce pressure-induced pain

Pulsation Therapy

- May help aid in the management of edema and patient comfort
- Helps improve capillary blood flow
- Aids in encouraging lymph flow which may reduce edema

Scheduled Pulsation

The system provides the convenience and safety of preset automatic pulsation. A complete pulsation cycle can be adjusted to 8, 16 or 32 minutes.

Pressure Relief Therapy

The system provides pressure relief therapy, which significantly reduces interface pressures when compared to a standard hospital bed mattress.

Reduced Heel Area Pressure

The cushions at the foot of the bed operate at a lower pressure level, which reduces heel area interface pressures. The cushion corresponding to the heel of the patient may also be removed to further reduce heel area pressure.

InstaFlate™ Function

The system is programmed with an InstaFlate function. When this function is activated through the control screen interface, the mattress will rapidly inflate to maximum firmness. This allows easy patient transfer and repositioning.

Touch Screen Interface

The *TheraKair Visio* MRS is controlled via a touch screen interface that gives details of system status and allows for easy and precise adjustment.

Individualized Patient Support

The *TheraKair Visio* MRS is made up of fifteen air cushions that run over the width of the mattress. The pressures within the cushions are automatically controlled at the therapy unit and can be adjusted to accommodate the weight and body shape of each patient, facilitating pulsation therapy and pressure relief therapy.

Quick Release Air Supply Connector

The air connection between the mattress and the therapy unit can be quickly and easily disconnected in the event that CPR (cardiopulmonary resuscitation) is required for the patient. When disconnected, the mattress will rapidly deflate and the therapy unit will go into standby mode. It will take approximately 10 seconds for the mattress to deflate sufficiently for CPR to commence.

Warmer

A warmer is incorporated into the therapy unit to allow adjustment of air temperatures entering the mattress and to help increase patient comfort.

Reusability

The *TheraKair Visio* MRS is reusable after proper care and cleaning. Minor leaks and holes caused by needle and pin pricks will not prevent reuse or cause deflation.

Indications

TheraKair Visio MRS is indicated for:

- Edema management
- Prevention or treatment of up to stage IV pressure ulcers
- Patients requiring a surface to assist pain management

Contraindications

Patient conditions for which the application of pressure relief therapy on *TheraKair Visio* MRS is contraindicated include:

- Spinal cord injury
- Cervical traction

Intended Care Setting

Intended care settings of the *TheraKair Visio* MRS are:

- Acute Care
- Post Acute Care

Connecting the System to Other Devices

Devices That May Be Attached to This System

None.

Devices to Which This System May Be Attached

The *TheraKair Visio* MRS can be fitted to most standard hospital beds.

The *TheraKair Visio* Therapy Unit can be fitted to the footboard of most hospital beds.

Risks and Precautions

Mattress Replacement

The *TheraKair Visio* MRS is a mattress replacement system and should never be fitted before first removing the existing mattress from the bed. Fitting the *TheraKair Visio* mattress on top of an existing mattress will raise the patient above the side rails of the bed, placing them at risk of falling out of the bed or becoming trapped between the mattress and the side rails.

Transfer

Standard precautions should be taken during patient transfer.

Side Rails and Restraints

WARNING: Use or non-use of restraints, including side rails, can be critical to patient safety.

Serious or fatal injury can result from the use (potential entrapment) or non-use (potential patient falls) of side rails or other restraints. See related **Safety Information**.

Patient Migration

Specialty surfaces have different shear and support characteristics than conventional surfaces and may increase the risk of patient movement, sinking and / or migration into hazardous positions of entrapment and / or inadvertent bed exit. Monitor patients frequently to guard against patient entrapment.

Use With Other Devices

All *TheraKair Visio* MRS components are designed to be used as single system. The *TheraKair Visio* Therapy Unit should only be used with the *TheraKair Visio* MRS. Any attempt to connect and use the therapy unit with any other brand or model of mattress or overlay may result in improper operation of equipment, possibly leading to increased risk of patient injury.

Protection Against Hazards

Fluids

Avoid spilling fluids on any part of the therapy unit.

If spills do occur:

- Disconnect the unit from the main electricity supply
- Clean fluids from the case

Ensure that there is no moisture in or near the power inlet, power switch and power plug before reconnecting the power supply.

- Check the operation of controls and other components in the area of the spill
- Perform applicable checkout procedures

Fluids remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and staff.

Prolonged contact with therapy unit may result in hazard to patient or caregiver.

Power Cable

The system should never be operated with a worn or damaged power cable. Should the power cable be found to be worn or damaged, contact your local Arjo subsidiary or agent for a replacement.

Safety Information

Patient Size and Weight

Maximum recommended patient weight is 135 kg

Minimum recommended patient weight is 23 kg

Maximum recommended patient height is 200 cm

Minimum recommended patient height is 137 cm

Air Intake

The air intakes for the *TheraKair Visio* Therapy Unit are at the rear and on the underside. To avoid the risk of blocking these intakes, it is recommended that the unit should only be mounted on the footboard of a bed frame. Should the bed not have a suitable footboard, it is recommended that the unit is placed on a solid stationary surface.

Side Rails / Patient Restraints

Whether and how to use side rails or restraints is a decision that should be based on each patient's needs and should be made by the patient and the patient's family, physician and caregivers, with facility protocols in mind. Caregivers should assess risks and benefits of side rail / restraint use (including entrapment and patient falls from bed) in conjunction with individual patient needs, and should discuss use or non-use with patient and / or family. Consider not only the clinical and other needs of the patient but also the risks of fatal or serious injury from falling out of bed and from patient entrapment in or around the side rails, restraints or other accessories. In the US, for a description of entrapment hazards, vulnerable patient profile and guidance to further reduce entrapment risks, refer to FDA's Hospital Bed System Dimensional and Assessment Guidance To Reduce Entrapment. Outside the US, consult the local Competent Authority or Government Agency for Medical Device Safety for specific local guidance. Consult a caregiver and carefully consider the use of bolsters, positioning aids or floor pads, especially with confused, restless or agitated patients. It is recommended that side rails (if used) be locked in the full upright position when the patient is unattended. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency. Monitor patients frequently to guard against patient entrapment.

When selecting a standard mattress, ensure the distance between top of side rails (if used) and top of mattress (without compression) is at least 8.66 in (220 mm) to help prevent inadvertent bed exit or falls. Consider individual patient size, position (relative to the top of the side rail) and patient condition in assessing fall risk.

Bed Frame

Always use a standard healthcare bed frame with safeguards or protocols that may be appropriate. It is recommended that bed and side rails (if used) comply relative to the mattress to help minimize any gaps that might entrap a patient's head or body. Inappropriate, incorrect, poorly fitted or poorly maintained bed rails and mattress replacements can potentially produce gaps that might entrap a patient's head or body.

Bed Height

To minimize the risks of falls or injury the patient surface should always be in the lowest practical position when the patient is unattended. Make sure area under and around bed frame is clear of objects, persons and parts of body before adjusting height.

Brakes

Caster brakes should always be locked once the bed is in position. Verify wheels are locked before any patient transfer to or from the bed.

Head of Bed Elevation

Keep head of bed as low as possible to help prevent patient migration.

Patient Entrance / Exit

Caregiver should always aid patient in exiting the bed. Make sure a capable patient knows how to get out of bed safely (and, if necessary, how to release the side rails) in case of fire or other emergency.

Skin Care

Monitor skin conditions regularly and consider adjunct or alternative therapies for high acuity patients. Give extra attention to any possible pressure points and locations where moisture or incontinence may occur or collect. Early intervention may be essential to preventing skin breakdown.

No Smoking in Bed

Smoking in bed can be dangerous. To avoid the risk of fire, smoking in bed should never be allowed.

Therapy Unit

- Only use a grounded power receptacle and the power cable supplied with the *TheraKair Visio* Therapy Unit. The power cable should be positioned to avoid a tripping hazard and / or damage to the cable. Arjo recommends placing the cable under the bed frame and attaching it to a mains wall socket by the head of the bed. The system should never be operated with a worn or damaged power cable. Should the power cable be found to be worn or damaged, contact Arjo or an Arjo authorized representative for a replacement.
- Only use replacement fuses that have the same specified rating. Using fuses with higher ratings could result in damage or injury.
- The therapy unit is a precision electronic product. Use care when handling or transporting. Dropping or other sudden impacts may result in damage to the therapy unit.
- Do not open the therapy unit – risk of electrical shock. Do not attempt to repair or service the therapy unit. Repairs and service should be conducted only by Arjo or an Arjo authorized representative. If the therapy unit is not functioning properly or has been damaged, unplug the unit and take it out of service immediately. Contact Arjo or an Arjo authorized representative for repair and service information.
- Do not place any objects or items, such as blankets, on or over the therapy unit.
- The air intakes for the *TheraKair Visio* Therapy Unit are on the back of the unit. To avoid the risk of blocking these intakes, it is recommended that the unit be mounted on the footboard of a bed frame. Should the bed not have a suitable footboard, it is recommended that the unit be placed on a solid stable surface.

Interference

Although the *TheraKair Visio* Therapy Unit conforms to the intent of IEC 60601-1-2 in relation to electromagnetic compatibility, all electrical equipment may produce interference. If interference is suspected, move equipment from sensitive devices or contact the manufacturer. (IEC 60601-1-2. Medical Electrical Equipment – Part 1: General Requirements for Safety, Amendment No. 2. Collateral Standard: Electromagnetic Compatibility Requirements and Tests). See the **Electromagnetic Emissions Information** section for detailed testing data.

General Protocols

Follow all applicable safety rules and institution protocols concerning patient and caregiver safety.

If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor. In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

End of Life Disposal

- Fabric material used on the mattresses or any other textiles, polymers or plastic materials etc. should be sorted as combustible waste.
- Mattresses at the end of life should be disposed of as waste according to the national or local requirements which may be landfill or combustion.
- Pump units have electrical and electronic components should be disassembled and recycled per Waste of Electrical and Electronic Equipment (WEEE) or in accordance with local or national regulation.

Customer Contact Information

For questions regarding this product, supplies, maintenance, or additional information about Arjo products and service, please contact Arjo or an Arjo authorized representative or see www.Arjo.com.

PREPARATION FOR USE

This chapter describes the steps necessary to prepare the *TheraKair Visio* MRS for use.

Check Out Procedures

The unit is shipped from Arjo ready to be put into service. Before putting the unit into service, the following steps should be performed.

1. Remove the mattress and the therapy unit from their packaging.

Do not use sharp instruments to open packaging. Damage could result to mattress, cover sheet, therapy unit or mattress hose.

2. Inventory and visually inspect that the following items are present and in good condition:
 - One *TheraKair Visio* MRS Therapy Unit
 - One *TheraKair Visio* MRS Form Fitted Cover Sheet
 - One *TheraKair Visio* MRS Mains Power Supply Cord
 - One *TheraKair Visio* MRS Mattress
 - One *TheraKair Visio* MRS Tote Bag
 - One *TheraKair Visio* MRS User Manual

Contact carrier and Arjo immediately if damage is observed.

Preparing the System for Use

Fitting the Mattress to the Bed

1. Level the bed and adjust to the lowest height.
2. Remove the hospital mattress from the bed.
3. Place the *TheraKair Visio* MRS Tote Bag on the hospital bed.
4. Open the tote bag and remove the mattress and the cover sheet. Unfasten the straps on the mattress. Place such that the air cushions are on top of the mattress and the hose connections are towards the foot of the bed.
5. Fit the cover sheet to the mattress.

Fitting the Mattress Cover Sheet

Fit the two head end corners of the cover sheet to the mattress.

Fit the closed foot end corner to the mattress.

Open the flaps on the remaining corner and wrap around the air supply hose.

Underside views


Attach the two hook and loop fasteners together to secure the corner.

Preparing the Therapy Unit

The therapy unit should be fitted to the footboard of the hospital bed, or if required, placed on a solid stationary surface.

1. Fit the therapy unit to a footboard of a hospital bed:
 - Pull out the green hanger at the back of the unit.
 - Fit over the footboard. Spring pressure will hold the unit in place.
2. Plug the mains power supply cord into the *TheraKair Visio* MRS Therapy Unit.

3. Plug power supply cord into a power outlet.

4. Plug the air supply connector into the therapy unit.

- i** The air supply connector will only connect to the therapy unit in one orientation. Rotate the connector until it locks in place.



OPERATION

This chapter describes how to correctly operate the *TheraKair Visio* MRS, how to operate and navigate the touch screen and how to safely power up and power down the system.

Power-Up Procedures

Initial Power-Up

1. Connect the system as described in the Preparation for Use section of this manual.
2. Press the green power switch located at the side of the therapy unit to the ON (I) position.

The indicator light on the top of the therapy unit should light up green.

Power-Up After an Improper Shutdown

Should the system be powered up following an improper shutdown of between 30 seconds and 30 minutes (power outage etc.), the system will start up and operate at the last stored settings for inflation and pulsation. Warming and InstaFlate function will be reset to zero.

If the shutdown is less than 30 seconds, the system will resume therapy at the last stored settings. Following any period of shutdown, the screen lock function will resume at the previous setting.

Power-Down Procedure

1. Press the green power switch located at the side of the therapy unit until it is in the OFF (O) position.
2. Unplug the therapy unit.

Performing CPR (Cardiopulmonary Resuscitation)

1. Disconnect the air supply connector by grasping the grip ring and pulling it away from the air supply unit.

It will take approximately 10 seconds for the mattress to deflate sufficiently for CPR to commence.

2. Swing the air supply hose away from the therapy unit.
3. Level the bed.
4. Begin CPR.

When the air supply connector is disconnected a CPR alarm screen will be displayed and the unit will go into standby mode.

After CPR has been performed and the patient is stable

1. Raise the side rails and lock into position.
2. Reconnect the air supply connector to the therapy unit.
3. Reset the alarm condition to resume therapy.

Control Screen Display

The control screen display on the front of the unit is a touch sensitive screen. The system control screens will be displayed on this screen. These screens will display information on current system operations and settings and allow the alteration of these values.

The screens have soft buttons; this means that the buttons operate exactly as real buttons would, simply press the buttons to activate them.

The operation of the control screens is detailed in the following pages.

The control screen should only be operated by finger; using pens or pointing devices will damage the screen and affect future use.

Control Screens

Home Screen

The following control screens can be accessed from this screen:

- Utilities menu screen
- Pulsation control screen
- Warming control screen
- InstaFlate function control screen
- Firmness control screen

There is a readout of current system settings at the bottom of each screen.

Pressing  on this or any subsequent screen will call up the quick reference guide providing help in operating the system.


Pressing and holding  for three seconds will lock / unlock the control screen display.

This button appears only on the home screen.

Pressing  on any subsequent screen will return to the previous screen.

Pressing  on any subsequent screen will confirm the settings entered on that screen and return to the previous screen.

Utilities Screen

The following control screens can be accessed from this screen:

- About screen
- Settings screen
- Alarms screen
- Language screen

About Screen


The About screen gives details of the version of software that the *TheraKair Visio* MRS is running.

Pressing and holding  for three seconds will call up the service screens via an access code screen.

This button appears only on the About screen.


Settings Screen

The Settings screen allows you to adjust the backlight intensity on the control screen display and to reset both the pulsation and therapy counters.

- Press the up or down arrows to change the screen brightness.
- Press and hold each reset button to reset and restart the relevant counter.


Alarms Screen

From the Alarms screen the volume of alarm warnings can be set.

- Press up and down to select the volume.
- Press test to sound the alarm at the set level.
- Press save to save current settings.
- Press reset to return alarm settings to the previously saved setting.


Language Screen

The language displayed on the control screens is selected on this screen.


- Press the up or down arrows to scroll through available languages.
- Press the save button to change to the selected language.
- Press reset to return the language settings to the previously saved setting.

Access Code Screen

The Access Code screen prevents unauthorized access to the service screens. The correct access code must be entered before proceeding.

Pulsation Screen


The Pulsation screen is accessed by pressing the pulsation button on the Home screen. This screen allows the setting of:

- The intensity of pulsation therapy to low, medium or high.
- The cycle time of the pulsation therapy to 8, 16 or 32 minutes

Pulsation therapy can also be toggled ON / OFF as required.

Warming Screen


The Warming screen allows the setting of the patient warming to low, medium or high. Warming can also be toggled ON / OFF as required.

InstaFlate Screen

- When the InstaFlate function is activated, the mattress will be rapidly inflated to its maximum firmness.
- From this screen the InstaFlate function can be switched ON or OFF.
- The exit button will only appear on this screen when the InstaFlate function is deactivated.
- The InstaFlate function will deactivate automatically after 15 minutes.
- When the InstaFlate function is cancelled, the system will resume therapy with the previous settings.
- Excessive use of the InstaFlate function is not recommended as this reduces patient comfort and may negate therapy.

Firmness Screen


The Firmness screen allows the setting of mattress firmness. Pressing the presets button will call up the patient presets screen where the patients height and weight can be entered.

Presets Screen

The patients height and weight are entered on the Presets screen. The *TheraKair Visio* MRS uses this information to estimate the level of firmness required. This screen can be accessed from the Firmness screen by pressing the Presets button.

Quick Reference Guide

The quick reference guide provides help in operating the system. This screen is accessed by pressing  on any other screen.

Alarms and Alerts

During operation, if a fault or situation that impairs functionality or therapy is detected by the system, an alarm will sound and an alarm message screen will be displayed. The status indicator light above the control screen will change to blue until the alarm is resolved.


The alarm message screen will give details of the particular alarm trigger and suggest a course of action.

When the cause of the alarm has been removed, press the reset button to resume therapy.

There are five alarm screens and two alert screens; they are as follows.

All alarms are classified as low priority alarms.


CPR Alarm Screen

This screen is displayed when the air supply connector is disconnected.

When the connector is reconnected, press the reset button to clear this alarm.

The alarm sounder may be paused for two minutes by pressing the pause button.

System Integrity Alarm

This screen is displayed when the system experiences a fatal error.

Contact your local Arjo subsidiary or agent for more details.

The alarm sounder may be paused for two minutes by pressing the pause button.

Air Pressure Alarm

This screen is displayed when the system air pressure falls outside the set parameters.

This may be due to damage to the mattress or to the therapy unit.

Check the mattress and hose for damage.

Contact your local Arjo subsidiary or agent for more details.

The alarm sounder may be paused for two minutes by pressing the pause button.

Alarm

Unit is over-temperature

Over Temperature Alarm

This screen is displayed when the therapy unit has overheated, this may be due to a blockage or a high ambient temperature.

Check that the air vents on the system are clear from obstruction and that the air filter is clean.

Power the system down and allow it to cool sufficiently.

Alarm

Therapy in-active

Therapy Inactive Alarm


This screen is displayed when the system has been left in an alert condition for more than two minutes.

The alarm sounder may be paused for two minutes by pressing the pause button.

Alert

Mattress not connected!

Connect mattress?

Mattress Not Connected Alert

This screen is displayed when the system is powered up before connecting the air supply connector to the therapy unit.

Connect the air supply connector and press the YES button to clear this alert.

Press NO to display the CPR Alarm screen.

Alert

Mattress not connected!

Ready to start?

This screen is displayed after pressing the YES button on the previous screen.

Connect the air supply connector and press YES to begin therapy and clear this alert.

PATIENT PLACEMENT

This chapter gives details of procedures for moving the patient safely to and from the *TheraKair Visio* MRS


Before placing a patient on the *TheraKair Visio* MRS, it is recommended that personnel familiarize themselves with the information outlined in the **Safety Information** section and **Operation** chapter of this manual as well as the procedures outlined below.

Preparing for Patient Placement

1. Ensure that the system is set up as described in the **Preparation for Use** chapter of this manual.
2. Ensure that a cover sheet is fitted and secured before transferring a patient.
3. Ensure that the system has been properly cleaned and disinfected, as detailed in the **Care and Cleaning** section of this manual before a new patient is moved to the *TheraKair Visio* MRS.
4. A clean set of air cushions and a clean cover sheet should be used for each new patient.

Patient Transfer to *TheraKair Visio* MRS

1. Press the InstaFlate button on the Home screen to access the InstaFlate screen.
2. Press the ON button on the InstaFlate screen to activate the InstaFlate function.
3. Raise / lower the height of the bed to the same level as the surface from which the patient is being transferred.
4. Lower the side rails.
5. Transfer the patient to the *TheraKair Visio* MRS.
6. Center the patient on the mattress from side to side.
7. Raise the side rails on both sides of the bed and lock in place.


8. Press the OFF button on the InstaFlate screen to deactivate the InstaFlate function.
9. Exit the InstaFlate screen to return to the Home screen.
10. Adjust the bed height to its lowest level and adjust the head of the bed as needed for patient comfort.
11. Make patient specific adjustments as required.

Patient Specific Adjustments


Foot End Interface

If necessary, heel area pressure can be reduced by disconnecting a foot end cushion corresponding to the patients heel. It is not recommended that more than one cushion be disconnected at a time.

1. Locate the cushion corresponding to the patients heel.
2. The cushions are connected to the air supply at the side opposite the air supply hose.

3. Disconnect the cushion by depressing the grey latch and pulling the connector out.


4. Cushions are reconnected by pushing the connector back into the socket.

Pulsation Adjustment

If necessary, pulsation therapy may be activated.

1. Press pulsation on the Home screen to call up the Pulsation screen.
2. Select the required pulsation intensity.
3. Select the required pulsation cycle time.
4. Perform a hand check to verify sufficient support:
 - Slide your hand flat between the mattress foam base and the cushions to verify a 2.5 cm to 4 cm clearance.


Warming Adjustment

If necessary warming may be used to increase patient comfort.

1. Press warming on the Home screen to call up the Warming screen.
2. Select the level of warming required.
3. Verify patient comfort.

Firmness Adjustment

1. Press the firmness button on the Home screen to enter the Firmness screen.
2. The firmness is adjusted by pressing the arrows up or down. With a new patient the presets must first be set as follows:
 - Press the presets button on the Firmness screen
 - Enter the patients weight and height
 - Press the set button to save the details and return to the Firmness screen
3. Perform a hand check to verify sufficient support:
 - Slide your hand flat between the mattress foam base and the cushions to verify a 2.5 cm to 4 cm clearance.

4. Allow sufficient time for the air supply unit and mattress to respond to air pressure changes.

When adjusting the patients position (e.g. from Fowler's position to supine) the firmness setting may need to be readjusted).

Patient Transfer from *TheraKair Visio* MRS

1. Press the InstaFlate button on the Home screen to access the InstaFlate screen.
2. Press the ON button on the InstaFlate screen to activate the InstaFlate function.
3. Adjust the bed height as required.
4. Lower the side rails.
5. Position a wheelchair or gurney next to the bed.
6. Transfer the patient.
7. Exit the InstaFlate screen to return to the Home screen.

Patient Transport on the *TheraKair Visio* MRS

Before transporting a patient the mattress should be fully deflated and the side rails raised. It is recommended that the period of time that the patient is resting on a deflated mattress be as short as possible.

NURSING CARE

This chapter gives details of procedures for general nursing care.

Bathing the Patient


1. To facilitate bathing, press the InstaFlate button on the Home screen to access the InstaFlate screen.
2. Press the ON button on the InstaFlate screen to activate the InstaFlate function.
3. Bathe patient.
4. Wipe cover sheet clean while bathing patient.

After Bathing the Patient

1. Center the patient on the *TheraKair Visio* MRS
2. If the InstaFlate function is still active, press the OFF button on the InstaFlate screen and exit the screen to resume therapy at previous settings.

Skin Care

1. Keep the skin dry and clean.
2. Check patient skin regularly, particularly in areas where incontinence and drainage occur.


Incontinence / Drainage

1. The *TheraKair Visio* MRS Cover Sheet should always be placed over the mattress.
2. Breathable underpads are recommended for incontinent patients as plastic backed underpads will block moisture vapor transmission and airflow from the mattress.
3. Watch for incontinence or drainage and provide meticulous skin care following each episode.

Mattress Firmness Adjustment

When adjusting the patients resting position (e.g. from Fowlers position to supine) it may be necessary to adjust the firmness setting.

1. Press the firmness button on the Home screen to enter the Firmness screen.
2. The firmness is adjusted by pressing the arrows up or down.
3. Perform a hand check to verify sufficient support.
4. Slide your hand flat between the mattress foam base and the cushions to verify a 2.5 cm to 4 cm clearance.

Bedpan Replacement

1. Press the pulsation button on the Home screen to enter the Pulsation Control screen.
2. Press the ON / OFF button to switch pulsation off.
3. Lower the side rail.
4. Position the bedpan parallel to the patients buttocks.
5. Turn the patient towards the opposite side rail.
6. Position the bedpan underneath the patient by pushing the bedpan down into the mattress.
7. Gently roll the patient onto the bedpan keeping one hand on the bedpan.
8. Raise the side rail.
9. Raise the head of the bed to a comfortable height for the patient.

Bedpan Removal

1. Lower the head of the bed to a flat position.
2. Lower the closest side rail.
3. Grasp the bedpan firmly with one hand and gently roll the patient to the opposite side.
4. Grasp the bedpan with both hands, push down into the mattress and slide out.
5. Center the patient on the mattress.
6. Raise the side rail.
7. Press the ON / OFF button to switch pulsation on.

CARE AND CLEANING

This chapter describes general care and cleaning of the *TheraKair Visio* MRS. Both scheduled and unscheduled maintenance are described along with preventative maintenance. Information on transport and storage of the system while not in use is also detailed.

The following are the Arjo recommended cleaning procedures for the *TheraKair Visio* MRS and periodic maintenance for the therapy unit while in use.

Daily Care and Cleaning of *TheraKair Visio* MRS While in Use

The patient need not be removed from the *TheraKair Visio* MRS when performing daily cleaning procedures.

Mattress

1. Press the InstaFlate function button on the Home screen to access the InstaFlate function screen.
2. Press the ON button on the InstaFlate function screen to activate InstaFlate function.
3. Verify the side rails on both sides of the bed are locked in the full upright position.
4. Gently roll the patient to one side of the mattress.
5. Wipe the cover sheet as required with soap and water or approved hospital germicide, according to individual hospital policy.
6. Extensively soiled cover sheet or mattress should be replaced and laundered. Refer to the **Weekly Care and Cleaning** section of this chapter for cleaning procedures.
7. Return the patient to the center of the mattress.
8. If InstaFlate function is still active, press the OFF button on the InstaFlate function screen and exit the screen to resume therapy at previous settings.

Therapy Unit Screen

The Therapy Unit screen should be cleaned daily, or more often as required.

The Control screen display should be locked before cleaning.

Pressing  on the Home screen, and holding for three seconds will lock / unlock the control screen display.

The screen may be wiped with a cloth dampened in water and mild detergent.

Use only non-abrasive cleaning products on the therapy unit and control screen.

Weekly Care and Cleaning of *TheraKair Visio* MRS While in Use

In addition to the daily care and cleaning, the therapy unit may be wiped with a cloth dampened in water and mild detergent. If required, the therapy unit may then be disinfected using a cloth dampened in chemicals recommended by your local Arjo subsidiary or agent.

TheraKair Visio MRS Cover Sheet

The *TheraKair Visio* MRS Cover Sheet should be replaced weekly, or as determined by caregiver. Excessively soiled mattresses and cover sheets should be replaced and laundered. Consult your local Arjo subsidiary or agent for further information on recommended chemicals and washing instructions.

Inflatable Air Cushions

Remove the *TheraKair Visio* MRS Cover Sheet and detach the air cushions from their fasteners at both sides of the mattress. Detach the air cushions from the air connectors and withdraw each air cushion from the mattress. Wipe each air cushion with water and mild detergent. If required, the air cushion may then be disinfected using chemicals recommended by your local Arjo subsidiary or agent.

Power Cable

The *TheraKair Visio* MRS should not be operated with a worn or damaged power cable. Inspect the power cable for any signs of wear or damage and contact your local Arjo subsidiary or agent for a replacement.

Periodic Maintenance of the Therapy Unit

Periodic maintenance to be performed at regular intervals while the unit is in operation and after each patient includes the following:

1. Operate the unit's functions and observe the operation of the unit as outlined in the **Operation** chapter of this manual. If any abnormalities are observed in the operation of the unit contact Arjo.
2. Inspect the mattress, cover sheet, therapy unit and air supply hose for wear and damage. Unserviceable items should be replaced with an identical item.
3. Check the air filter, normally every week while the unit is in service. The filter is in a drawer that opens to the front of the case:
 - Power down the system.
 - Unlatch the drawer by lifting the button under the 'handle' part of the drawer.
 - Slide the drawer fully out.
 - Raise the front flaps on the drawer to access the filter.
4. If required, the filter may either be replaced or washed with a mild detergent and allowed to dry naturally before refitting.

Safe Transport and Storage of the *TheraKair Visio* MRS

Transport and Storage of Individual System Parts

It is recommended that the system should be transported and stored in the packaging provided. It is also recommended that the *TheraKair Visio* MRS be stored in the bag provided.

Transport of System While Fitted to Hospital Bed

When transporting the system while attached to a hospital bed, ensure that the air supply hose and the power cord are safely away from the wheels of the bed.

Take care that the therapy unit does not get damaged while pushing the bed through hospital doors.

SPECIFICATIONS

Specifications subject to change without notice.

This chapter gives the various physical and technical specifications of the *TheraKair Visio* MRS and details of spare and replacement parts.

Maximum Recommended Patient Weight	135 kg
Minimum Recommended Patient Weight	23 kg
Maximum Recommended Patient Height	200 cm
Minimum Recommended Patient Height	137 cm

Classification

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

Continuous Operation

- Type B Equipment
- Class I Equipment
- Ordinary Equipment

Mattress Overlay

Dimensions:

(Inflated) (90 cm)	208 cm x 90 cm x 20.3 cm
(Inflated) (80 cm)	208 cm x 80 cm x 20.3 cm
Foam Pad (90 cm)	208 cm x 90 cm x 7.6 cm
Foam Pad (80 cm)	208 cm x 80 cm x 7.6 cm
Weight	13.2 kg

Therapy Unit

Dimensions	41.2 cm x 26.4 cm x 21.0 cm
Weight	6.3 kg

Electrical Data

Voltage	100 - 240 VAC
Frequency	50 Hz - 60 Hz
Current	3.7A (100V) / 2.3A (240V)
Fuse rating	1500A IC, 250V 6.3A
Fuse Type	Time lag (Slo-Blo)

Electrical Data - Kingdom of Saudi Arabia

Voltage	230 VAC
Frequency	60 Hz
Current	2.3A
Fuse rating	1500A IC, 250V 6.3A
Fuse Type	Time lag (Slo-Blo)

Transport and Storage Conditions

Storage Temperature Range	-10°C to +50°C
Relative Humidity Range	10% to 85% including condensation
Barometric Pressure Range	700 hPa to 1060 hPa

Operating Conditions

Operational Temperature Range	+10°C to +35°C
-------------------------------------	----------------

The life expectancy of the cover sheet and cushions is up to 25 cycles (depending on washing and cleaning procedures). The TheraKair Visio MRS conforms to the following International Standards: EN60601-1:2006/A1:2013; ANSI/AAMI ES60601-1(2005)+AMD(2012) and CAN/CSA-C22.2 No.60601-1(2008)+(2014)

The *TheraKair Visio* MRS is classified as a Type B applied part under IEC 60601-1:2005/A1:2012.

REPLACEMENT PART NUMBERS

<i>TheraKair Visio</i> MRS Therapy Unit	M6257763
<i>TheraKair Visio</i> MRS Therapy Unit	413412
<i>TheraKair Visio</i> MRS Therapy Unit (Kingdom of Saudi Arabia)	415165
<i>TheraKair Visio</i> MRS User Manual	340750-AH
<i>TheraKair Visio</i> MRS Power Cord	
UK/KSA	M4268743
Europe	M4268742
USA, Japan, Canada	M4268840
Australia / New Zealand	M4268893
<i>TheraKair Visio</i> MRS Mattress:	
800 mm wide	M6265529
900 mm wide	M6265526
<i>TheraKair Visio</i> MRS Tote Bag:	
900 mm wide	M6248435
<i>TheraKair Visio</i> MRS Air Filter	
Small Filter	M3248823
Large Filter	M3248824
<i>TheraKair Visio</i> MRS Cover Sheets	
800 mm wide	M6246643
900 mm wide	M6246635
<i>TheraKair Visio</i> MRS Air Cushions	
800 mm wide	
Head End	M6240535
Foot End	M6240536
900 mm wide	
Head End	M6240531
Foot End	M6240532
<i>TheraKair Visio</i> MRS Mattress Foam	
800 mm wide	M3249386
900 mm wide	M3249376
<i>TheraKair Visio</i> MRS Therapy Unit Case	M3275033

ELECTROMAGNETIC COMPATIBILITY (EMC)

Product has been tested for compliance with current regulatory standards regarding its capacity to block EMI (electromagnetic interference) from external sources.

Some procedures can help reduce electromagnetic interferences:

- Use only Arjo cables and spare parts to avoid increased emissions or decreased immunity which can compromise the correct functioning of the equipment.
- Ensure that other devices in patient-monitoring and/or life-support areas comply to accepted emissions standards.

WARNING: Wireless communications equipment such as wireless computer network devices, mobile phones, cordless telephones and their base stations, walkie-talkies, etc. can affect this equipment and should be kept at least 1.5m away from the equipment.

Intended Environment: Home Healthcare Environment and Professional Healthcare facility environment.

Exceptions: HF Surgical Equipment and the RF Shielded room of an ME SYSTEM for magnetic resonance imaging

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Guidance and manufacturer’s declaration - electromagnetic emission		
Emissions test	Compliance	Guidance
RF emissions CISPR 11	Group 1	This equipment uses RF energy only for its internal functions. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	This equipment is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electromagnetic Discharge (ESD) IEC 61000-4-2	±2kV, ±4kV, ±8kV, ±15kV air ±8kV contact	±2kV, ±4kV, ±8kV, ±15kV air ±8kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity level should be at least 30%
Conducted disturbances induced by RF fields EN 61000-4-6	3V in 0,15 MHz to 80 MHz 6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	3V in 0,15 MHz to 80 MHz 6V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than 1.0m, if the transmitter's output power rating exceeds 1W ^a Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with this symbol:
Radiated RF electromagnetic field EN 61000-4-3	Home Healthcare environment 10 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz	Home Healthcare environment 10 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz	
Electrical fast transient/burst EN 61000-4-4	±1kV SIP/SOP ports ±2kV AC port 100 kHz repetition frequency	±1kV SIP/SOP ports ±2kV AC port 100kHz repetition frequency	Mains power supply should be that of a typical commercial or hospital environment.
Power frequency Magnetic field EN 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity			
Surge	$\pm 0.5\text{kV}, \pm 1\text{kV}, \pm 2\text{kV}$, AC Mains Line to Ground	$\pm 0.5\text{kV}, \pm 1\text{kV}, \pm 2\text{kV}$, AC Mains Line to Ground	Mains power quality should be that of a typical commercial or hospital environment.
IEC61000-4-5	$\pm 0.5\text{kV}, \pm 1\text{kV}$, AC Mains Line to Line	$\pm 0.5\text{kV}, \pm 1\text{kV}$, AC Mains Line to Line	
Voltage dips, short interruptions and voltage variations on power supply input lines	0% U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase: at 0° 0% U_T ; 250/300 cycle	0% U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 cycles Single phase: at 0° 0% U_T ; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the pump requires continued operation during power mains interruption, it is recommended that the pump be powered from an uninterruptible power supply or battery.
NOTE: U_T is the AC mains voltage prior to application of the test level.			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1 V/m.</p>			

DEFINITION OF SYMBOLS USED

Warning of possible hazard to system, patient or staff

Important operational information

Tumble dry at low heat

May be dry-cleaned with solvent classified with the letter P

Do not bleach

Alternating current

On

Separate electrical and electronic components for recycling in accordance with the European Directive 2012/19/EC (WEEE).

Manufacturer name and address.

ETL Listed, Conforms to AAMI STD. ES60601-1, IEC STD 60601-1 Certified to CSA STD. C22.2 No.60601-1

Temperature Limitations

Humidity Limitations

Type B Applied Part

Do not iron

Wash temperature

Consult Instructions for Use

Off

2797

CE marking indicating conformity with European Community harmonised legislation. Figures indicate Notified Body supervision.

Read the IFU before use.

Indicates the product is a Medical Device according to EU Medical Device Regulation 2017/745

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AUSTRALIA
Arjo Australia Pty Ltd
78, Forsyth Street
O'Connor
AU-6163 Western Australia
Tel: +61 89337 4111
Free: +1 800 072 040
Fax: +61 89337 9077

BELGIQUE / BELGIË
Arjo NV/SA
Evenbroekveld 16
BE-9420 ERPE-MERE
Tél/Tel: +32 (0) 53 60 73 80
Fax: +32 (0) 53 60 73 81
E-mail: info.belgium@arjo.be

BRASIL
Arjo Brasil Equipamentos Médicos Ltda
Rua Marina Ciufui Zanelice, 329 PB02
Galpão - Lapa
São Paulo - SP - Brasil
CEP: 05040-000
Phone: 55-11-3588-5088
E-mail: vendas.latam@arjo.com
E-mail: servicios.latam@arjo.com

CANADA
Arjo Canada Inc.
90 Matheson Boulevard West
Suite 300
CA-MISSISSAUGA, ON, L5R 3R3
Tel/Tél: +1 905 238 7880
Free: +1 800 665 4831 Institutional
Free: +1 800 868 0441 Home Care
Fax: +1 905 238 7881
E-mail: info.canada@arjo.com

ČESKÁ REPUBLIKA
Arjo Czech Republic s.r.o.
Na Strži 1702/65
140 00 Praha
Czech Republic
Phone No: +420225092307
e-mail: info.cz@arjo.com

DANMARK
Arjo A/S
Vassingerødvej 52
DK-3540 LYNGE
Tel: +45 49 13 84 86
Fax: +45 49 13 84 87
E-mail:
dk_kundeservice@arjo.com

DEUTSCHLAND
Arjo GmbH
Peter-Sander-Strasse 10
DE-55252 MAINZ-KASTEL
Tel: +49 (0) 6134 186 0
Fax: +49 (0) 6134 186 160
E-mail: info-de@arjo.com

ESPAÑA
Arjo Ibérica S.L.
Parque Empresarial Rivas Futura, C/
Marie Curie 5
Edificio Alfa Planta 6 oficina 6.1.-62
ES-28521 Rivas Vacia, MADRID
Tel: +34 93 583 11 20
Fax: +34 93 583 11 22
E-mail: info.es@arjo.com
FRANCE
Arjo SAS

2 Avenue Alcide de Gasperi
CS 70133
FR-59436 RONCQ CEDEX
Tél: +33 (0) 3 20 28 13 13
Fax: +33 (0) 3 20 28 13 14
E-mail: info.france@arjo.com

HONG KONG
Arjo Hong Kong Limited
Room 411-414, 4/F,
Manhattan Centre,
8 Kwai Cheong Road,
Kwai Chung, N.T.,
HONG KONG
Tel: +852 2960 7600
Fax: +852 2960 1711

ITALIA
Arjo Italia S.p.A.
Via Giacomo Peroni 400-402
IT-00131 ROMA
Tel: +39 (0) 6 87426211
Fax: +39 (0) 6 87426222
E-mail: Italy.promo@arjo.com

MIDDLE EAST
Arjo Middle East FZ-LLC
Office 908, 9th Floor,
HQ Building, North Tower,
Dubai Science Park,
Al Barsha South
P.O Box 11488, Dubai,
United Arab Emirates
Direct +971 487 48053
Fax +971 487 48072
Email: Info.ME@arjo.com

NEDERLAND
Arjo BV
Biezenwei 21
4004 MB TIEL
Postbus 6116
4000 HC TIEL
Tel: +31 (0) 344 64 08 00
Fax: +31 (0) 344 64 08 85
E-mail: info.nl@arjo.com

NEW ZEALAND
Arjo Ltd
34 Vestey Drive
Mount Wellington
NZ-AUCKLAND 1060
Tel: +64 (0) 9 573 5344
Free Call: 0800 000 151
Fax: +64 (0) 9 573 5384
E-mail: nz.info@Arjo.com

NORGE
Arjo Norway AS
Olaf Helsetts vei 5
N-0694 OSLO
Tel: +47 22 08 00 50
Faks: +47 22 08 00 51
E-mail: no.kundeservice@arjo.com

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Address page - REV 24: 04/2019

At Arjo, we are committed to improving the everyday lives of people affected by reduced mobility and age-related health challenges. With products and solutions that ensure ergonomic patient handling, personal hygiene, disinfection, diagnostics, and the effective prevention of pressure ulcers and venous thromboembolism, we help professionals across care environments to continually raise the standard of safe and dignified care. Everything we do, we do with people in mind.

ArjoHuntleigh AB
Hans Michelsensgatan 10
211 20 Malmö, Sweden
www.arjo.com

arjo

340750-AH

CE
2797

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Arjo Inc.